Instructions for use

MK-dent Air Motors

Version: 2025-12-23-RevB

Language Index

English Name	Native Name	ISO 639-1	Page
English	English	en	3
French	Français	fr	23
Spanish	Español	es	43
German	Deutsch	de	63
Italian	Italiano	it	83
Chinese	中文	zh	103
Swedish	Svenska	sv	124
Hungarian	Magyar	hu	144
Slovak	Slovenčina	sk	164
Czech	Česky	cs	184
Turkish	Türkçe	tr	204
Norwegian	Norsk	no	224
Bulgarian	Български	bg	244
Lithuanian	Lietuvių	lt	264
Portuguese (Portugal)	Português	pt	284
Arabic	العربية	ar	304
Greek	Αγγλικά	el	325

Table of contents

1	About these instructions for use	5
1.1	Explanation of warnings	5
1.2	Compliance	5
2	Intended purpose.....	6
2.1	Intended use.....	6
2.2	Indications	6
2.3	Contraindications	6
2.4	Intended patient groups.....	6
2.5	Intended users.....	6
2.6	User environment.....	6
3	Safety and requirements	7
3.1	General safety regulations	7
3.2	Reporting obligations	7
3.3	Side effects	7
3.4	Risk patient groups.....	7
4	Product description	8
4.1	Description	8
4.2	Combination with other products	8
4.3	Product components	8
4.4	Scope of delivery	8
4.5	Required materials that are not included	8
4.6	Configurations	8
5	Accessories, consumables and spare parts.....	9
5.1	Accessories	9
5.2	Consumables	9
5.3	Spare parts	9
6	Installation.....	10
6.1	Requirements for the environment	10
6.2	Unpacking.....	10
6.3	Assembly	10
6.4	Preparing the connection	10
6.4.1	Connection to the treatment unit.....	11
6.4.2	Connection to the instruments.....	11
6.5	Startup and operation	12
6.5.1	Check for correct operation	12

6.5.2	Setting the compressed air	12
6.5.3	Speed control and change of direction of rotation	12
6.5.4	Connecting the external water supply (AM0012, AM0014)	12
6.5.5	Water regulation (AM1116)	12
7	Treatment	13
7.1	Patient, operator and assistant protection	13
7.2	Disconnecting after use	13
8	Processing according to ISO 17664-1	14
8.1	Pre-treatment at the point of use	14
8.2	Manual processing	15
8.2.1	Manual cleaning of the surface	15
8.2.2	Manual interior cleaning and drying	15
8.2.3	Manual disinfection	16
8.3	Automatic processing	16
8.4	Sterilization	16
8.5	Storage	17
9	Maintenance	17
9.1	Regular maintenance and checks	17
9.2	Cleaning in case of problems	17
10	Transport and disposal	18
11	Troubleshooting	18
12	Repair	19
12.1	Replacing the O-rings	19
12.2	Replacing the light source	19
12.2.1	Replacing the light source for AM1116	19
12.2.2	Replacing the light source for AM1006	19
13	Technical data	20
14	Explanation of symbols	20
15	Guarantee and warranty	21
15.1	Guarantee	21
15.2	Legal warranty	21
15.3	Conditions	21
15.4	Utilization	22
15.5	Exclusion of liability	22

1 About these instructions for use

These instructions for use apply to all MK-dent air motors (see section **4 Product description**).

All products listed here are only intended for dental treatment in the field of dentistry. Any kind of misuse or modification of the product is not permitted and may lead to hazards. All products listed here may only be used by qualified users in accordance with **2.5 Intended users** for the intended purpose and in compliance with the instructions for use.

1.1 Explanation of warnings

These instructions for use are intended to describe the safe and effective operation of the product to the user. Before operating the product, you must read these instructions for use, observe all safety instructions and warnings and strictly adhere to them. Pay particular attention to all information and procedures described in section **3 Safety and requirements**.

These instructions for use are part of the product and must be kept within reach of the product so that they are always accessible.

To avoid personal injury and damage to property, observe the warnings and safety information contained in this document. These are specially marked:

WARNING

A WARNING alerts you to a potentially serious outcome, unwanted result or safety risk. Failure to observe a warning can result in the death or serious injury of the user or patient.

CAUTION

CAUTION alerts you to situations requiring particular care for the safe and effective use of the product. Disregard for precautionary measures can result in minor or moderately severe personal injuries or damage to the product or other property and may potentially lead to a remote risk of serious injuries and/or environmental pollution.

NOTE

NOTE

This identifies special tips, such as help to the user or for the improvement of a work process.

1.2 Compliance

If you have any further questions regarding the applicable national or international standards, please contact the manufacturer:

MK-dent GmbH, Marie-Curie-Str. 2, 22941 Bargteheide,
Germany
+49 4532 40049-0
info@mk-dent.com

2 Intended purpose

2.1 Intended use

Air motors are an accessory for dental handpieces without direct patient contact, which are used in the preparation of dentures and for polishing teeth and restoration surfaces because they provide the necessary drive.

2.2 Indications

Not applicable. Air motors are an accessory for dental straight and contra-angle handpieces without direct patient contact.

2.3 Contraindications

There are currently no known contraindications.

2.4 Intended patient groups

Not applicable. Air motors are an accessory for dental straight and contra angle handpieces without direct patient contact.

2.5 Intended users

An air motor is only intended for use by licensed dentists or dental hygienists with appropriate training. Secondary users are the practice staff who are responsible for transport, processing and cleaning.

2.6 User environment

Application takes place in dental practices or dental facilities. Users wear protective gloves. The peripherals used fulfil the requirements of ISO 7494.

3 Safety and requirements

3.1 General safety regulations

WARNING

- This product is supplied non-sterile. Disinfect it before first use according to section **8 Reprocessing according to ISO 17664-1**.
- For proper reprocessing, reprocess the device after each patient in accordance with section **8 Reprocessing according to ISO 17664-1**. Safe and effective reprocessing after multiple uses without reprocessing the device leads to mixed infections and prevents proper reprocessing due to atypical accumulation of contaminants, occlusions and cleaning agent residues – especially in cavities and supply lines.
- This product may be used only by qualified personnel and only for dental handpieces.
- Never disassemble or modify this product.
- This instrument must not be repaired by the user in the event of a malfunction. Contact your dealer or the MK-dent repair service.
- Please do not use the product if the performance of the product deteriorates or in case of malfunctions.
- Do not use this product if it has not been serviced properly.
- Observe the precautionary measures for use in order to ensure correct and safe functioning.

CAUTION

- Before initial use and after each treatment, the instrument and accessories must be prepared in accordance with the instructions for use.
- The product may only be used in the field of dentistry with a dental treatment unit that has a CE marking or is otherwise certified/approved according to the valid national legal regulations and that is provided with a suitable supply of water and compressed air in compliance with the applicable standards.
- Test for proper functioning before each use and inspect the instrument visually for external damage and contamination.

3.2 Reporting obligations

In the event of serious deterioration in the patient's health during the use of this product, please report this to your specialized dealer or to us as the manufacturer – and to the national competent authority.

Our contact details are included in these operating instructions. The contact details of the notified body will also be provided to you upon request.

3.3 Side effects

There are currently no known side effects.

3.4 Risk patient groups

There are currently no known risk patient groups.

4 Product description

4.1 Description

Air motors are an accessory for dental handpieces that are used in the preparation of dentures and for polishing teeth and restoration surfaces because they provide the necessary drive.

Thanks to the quick-coupling system, replacing air-operated dental handpieces is easy, enabling immediate use. No complicated installation is required, allowing for simple integration in your treatment system.

Driven by compressed air and water from the treatment unit, an air motor does not affect other medical devices, nor is it affected by electromagnetic waves.

4.2 Combination with other products

REF MK-dent air motor	AM0002, AM0004, AM0012, AM0014, AM1002, AM1004, AM1006, AM1116
REF MK-dent straight and contra-angle handpiece	LB01, LB01W, LB02, LB1111, LB1111C, LB11W, LB4111, LB4111P, LE01, LE01L, LE11, LE11L, LE148, LE148L, LE148LP, LE148P, LE15, LE15L, LE16, LE16L, LE21, LE21L, LE21P, LE27, LE27L, LE54, LE54L, LE54LP, LE54P, LE74, LE74L, LC01, LC01L, LC11, LC11L, LC27, LC27L, LC74, LC74L, LP01L, LP11L, LP15L, LP27L, LP54L, LP74L
Connection*	ISO connection*

* The straight and contra angle handpieces are compatible with all motors with a connection in accordance with DIN EN ISO 3964.

4.3 Product components

Not applicable

4.4 Scope of delivery

The scope of delivery of the following components and accessories is included in the delivery unit:

REF MK-dent air motor	LT1014P Oil adapter	SP0010 External water connection
AM0012, AM0014	X	X
AM0002, AM0004, AM1002, AM1004, AM1006, AM1116	X	--

When you open the delivery package for the first time, make sure that the seal on the packaging is undamaged and that all components are included in the scope of delivery. If anything is missing, contact your dealer immediately.

4.5 Required materials that are not included

NOTE

NOTE

The scope of delivery does not include a dental handpiece.

The scope of delivery does not include MK-dent Premium Service Oil LU1011.

4.6 Configurations

Not applicable

5 Accessories, consumables and spare parts

5.1 Accessories

Not applicable

5.2 Consumables

Not applicable

5.3 Spare parts

CAUTION

The use of unauthorized accessories or unauthorized modifications to the product can lead to injuries or instrument malfunctions.

→ Only use spare parts and components that have been approved by MK-dent for combination with the product.

AM0002, AM0004	
Replacement motor without light	AM0002, AM0004
Replacement O-ring	2 x OR3021

AM0012, AM0014	
Replacement motor without light	AM0012, AM0014
Replacement O-ring	2 x OR3021

AM1002, AM1004, AM1006	
Replacement motor without light	AM1002, AM1004
Replacement motor with LED	AM1006
Replacement LED	BU8012M1
Replacement O-ring	3 x OR3021

AM1116	
Replacement motor with LED	AM1116
Replacement LED	BU8012KM
Replacement O-ring	3 x OR3011

6 Installation

6.1 Requirements for the environment

To supply the air motors with clean and dry air, clean the compressor air filter in your dental office. Do not supply moist compressed air. For drainage and regular maintenance of your compressor, please read the instructions for use for your air compressor. Ensure that the compressor is maintained in accordance with its instructions for use.

6.2 Unpacking

When opening the package for the first time, make sure that the content according to **4.4 Scope of delivery** is enclosed with the product. Make sure that the air motor has no defects or damage such as cracks.

6.3 Assembly

Not applicable

6.4 Preparing the connection

WARNING

When connecting your handpieces to this product, check that the connection is secure. Disconnecting the instrument from the motor during a treatment can be dangerous for patients and users!

Instruments must not be put on or taken off while the motor is rotating. This can cause damage to the motor and the instrument!

Never leave the motor running without cooling air, as this can lead to overheating.

Check for the absence of any damage, cracks, breakage or corrosion on the product and that it is safe to use.


CAUTION

The air motors must be disinfected before initial use and after each treatment. Only use dry, clean compressed air. Contaminated or moist compressed air leads to premature wear.

When a handpiece is removed, water and oil remain on the motor connection. Remove the water and oil completely. Make sure that no water or oil remains in your connection hose. If present, remove with a dry cloth.

6.4.1 Connection to the treatment unit

Attach the hose connection to the air motor and tighten the thread.

6.4.2 Connection to the instruments

Lightly oil the O-rings on the motor with MK-dent Premium Service Oil LU1011.

Push the instrument onto the air motor until it clicks into place.

To remove, pull the instrument off the motor in the axial direction.

6.5 Startup and operation

6.5.1 Check for correct operation

WARNING

When connecting your handpiece to this product, check that the connection is secure. An unstable connection will result in injury due to the body and other parts falling off.

Press the foot switch of your treatment unit once or twice for about one second and check – when air and water are sprayed – whether the water supply is sufficient and no leakage of air or water can be detected at the connection of the motor.

If you have problems, read section **11 Troubleshooting**. Use the water volume regulator on your treatment unit to set the water supply.


6.5.2 Setting the compressed air

WARNING


The compressed air at the motor must be set in accordance with the handpiece's instructions for use.

6.5.3 Speed control and change of direction of rotation

Operate the rotary control to control the speed. Turning the control knob continuously adjusts the speed in both clockwise and counter-clockwise rotation. If the markings on the rings are on top of each other, the motor will not rotate.

6.5.4 Connecting the external water supply (AM0012, AM0014)

Attach one end of the plastic hose to the water outlet valve of the air motor and the other end to the water inlet of the handpiece / contra-angle handpiece. Switch on the water supply to your treatment center. As the hose supplied is longer than necessary, cut it to the appropriate length before use.

6.5.5 Water regulation (AM1116)

The amount of water in the spray can be adjusted by turning the regulating ring at the rear of the air motor.

7 Treatment

7.1 Patient, operator and assistant protection

WARNING

Before the first contact with the patient, the product must be properly processed to prevent infection. See section **8 Processing according to ISO 17664-1**.

Protect the patient accordingly before use.

During treatment, the operator and assistants must wear gloves, protective goggles and masks to prevent infection.

CAUTION

Always stop working whenever operating noises become louder or irregular, there are stronger vibrations or in any case of overheating or visual damage. In these cases, contact your dealer or the MK-dent repair service.

7.2 Disconnecting after use

WARNING

The product must not be dropped.

Disconnect the air motor in accordance with your handpiece's instructions for use. When reconnecting, first clean the motor connection, connection hose and the inside of the connection in accordance with

6.4 Preparing the connection.

8 Processing according to ISO 17664-1

NOTE

NOTE

The possible number of applications does not depend on the processing cycles, but rather on the use of the device and the resulting wear and tear. The product's service life is therefore determined by signs of wear and tear and damage occurring during use.

Dispose of damaged products in accordance with section **10 Transport and disposal**. Do not use any products with obvious signs of wear or damage.

Only use cleaning agents / disinfectants that are authorized for the processing of medical devices. Follow the manufacturer's instructions for the cleaning agents / disinfectants.

If unsuitable cleaning agents or disinfectants are used, this can have negative effects on the instruments:

- Damage or corrosion
- Discoloration of the product
- Corrosion of metal parts
- Reduced service life

Follow the instructions of the cleaning agent / disinfectant manufacturer. Only use non-fixing cleaning agents and disinfectants.

Motors cannot be sterilized. Sterilization is not absolutely necessary for the intended use.

Ultrasonic cleaning is not recommended for processing the instruments, as this could reduce the product's service life.

8.1 Pre-treatment at the point of use

WARNING

Risk of infection from non-sterile instruments. Always wear protective gloves.

CAUTION

The use of disinfectant baths and disinfectants containing chlorine can lead to defects and malfunctions of the instrument.

The cleaning result is diminished by surface drying of residues.

→ Flush the instrument immediately when you finish working on the patient.

Do not place the instrument in solutions or similar.

NOTE

NOTE

Clean the connection hose and the inside of the instrument connector according to **6.4 Preparation for use** before reconnecting the instrument.

- Pre-cleaning of the instruments should be carried out immediately after treatment, taking personal protection into account. The aim is to prevent organic material and chemical residues from drying on the motor and to avoid contaminating the environment.
- Remove adhering contaminants with a disposable cloth / paper towel.
- Rinse the product with water after use (temperature below 35 °C (95 °F)).
- Only use non-fixing cleaning agents or water at a temperature below 35 °C (95 °F), as there may otherwise be increased adhesion of contaminants, making further processing steps more difficult.
- Clean the air motor within one hour after each treatment to prevent it from drying out.
- Ensure that the instrument is stored safely and transported in dry condition in a closed container to the processing site.

8.2 Manual processing

CAUTION

Insufficient cleaning of the surface and of water and air supply lines.

Thorough manual pre-cleaning must be carried out to ensure successful cleaning.

The use of manual pre-cleaning requires training of the employees who carry out this processing instruction. This is necessary to ensure successful cleaning.

Malfunctions can occur if cleaning is done in the ultrasonic machine

→ Clean the instrument exclusively as described below.

8.2.1 Manual cleaning of the surface

Tools:

- Drinking water¹ 30 °C ± 5 °C (77 °F - 95 °F)
- Brush or medium-hard toothbrush

Brush the air motor under running drinking water¹.

8.2.2 Manual interior cleaning and drying

Flush the spray air and spray water line with drinking water¹.

Blow the motor with medical compressed air from the outside and from the inside until no more drops of water are visible.

¹ Drinking water quality according to the EU Drinking Water Regulation (total bacterial count max. 100 CfU/ml)
MK-dent-Instructions-for-Use-Airmotors-Multilingual-RevB
Version: 2025-12-23-RevB

8.2.3 Manual disinfection

WARNING

Air motors cannot be sterilized.

The use of disinfection baths and disinfectants containing chlorine can lead to defects and malfunctions in the motor.

→ Only clean and disinfect the air motor manually.

NOTE

NOTE

Use alcohol-based instrument disinfectant with a short residence time that does not cause any relevant adhesion of proteins on the handpieces. The instrument disinfectant must be bactericidal, fungicidal and antiviral according to RKI requirements. Never use disinfectants containing chlorine, as these can damage the handpieces. Do not use aldehydic disinfectants, as they are incompatible with the material. Only use disinfectants whose microbiological efficacy is confirmed or guaranteed by the manufacturer (e.g., VAH/DGHM registration or CE marking).

Tools:

- Instrument disinfectant
- Disposable wipe

Spray the disinfectant onto a disposable cloth and then wipe the air motor. Do not place the instrument in disinfectant baths!

Observe the exposure time specified by the disinfectant manufacturer. Follow the disinfectant's instructions for use.

8.3 Automatic processing

Not applicable

WARNING

Automatic reprocessing is not to be used for air motors.

8.4 Sterilization

Not applicable

WARNING

Sterilization in an autoclave is not to be used for air motors.

8.5 Storage

After disinfection and drying, check that there are no drops of residual water. Store the processed instruments in a dry, dark and cool room, protected from dust and with the lowest possible microbial load.

Temperature	-20 °C - +70 °C (-4 °F - +158 °F)
Relative humidity	5 % - 95 %, non-condensing
Air pressure	700 hPa - 1060 hPa (10 psi - 15 psi)

CAUTION

Instrument malfunction during operation after heavily cooled storage conditions
→ Allow heavily cooled instruments to warm up to room temperature before starting to use them.

The instrument must be completely dried before storage. Any remaining water can result in discoloration or rust.

Do not store the instrument in a place where it could be exposed to water.

Do not store the instrument in places where chemicals are stored or where corrosive gases exist.

9 Maintenance

9.1 Regular maintenance and checks

WARNING


Inappropriate maintenance or service can lead to malfunctions and/or deterioration.

→ Always conduct proper maintenance and servicing.

CAUTION

For safe use of the product, perform checks and maintenance of the instrument according to these instructions for use.

Only use MK-dent Premium Service Oil LU1011 for maintenance.

Remove the motor from the hose and spray into the smaller of the two connection holes for 2 seconds (see illustration)

Use the LT1014P oil adapter for this.

The air motor should be maintained in this way at least after every disinfection.

9.2 Cleaning in case of problems

In case of problems, clean each part of the instrument. The instrument must be disconnected from the hose for cleaning. Instructions can be found in chapter **8.1 Pre-treatment at the point of use**.

10 Transport and disposal

CAUTION

Users or third parties may be infected by contaminated medical devices. Process the product before transport or disposal or dispose of it properly in compliance with the applicable regulations for the disposal of contaminated medical waste.

Reprocess the instrument in accordance with section **8 Processing in accordance with ISO 17664-1** and return it to us in its original packaging. We will take care of the professional disposal of the device and the original packaging for you.

11 Troubleshooting

CAUTION

This instrument is not intended to be repaired by the user in the event of a malfunction. Contact your dealer or the MK-dent repair service.

Process this instrument according to section **8 Processing according to ISO 17664-1** and send it to your dealer or directly to MK-dent if it is defective. Please contact MK-dent customer service to arrange for shipping or collection for factory repair. Questions about noise, maintenance, processing and functionality can often be clarified in a conversation, making it unnecessary to send in the product.

12 Repair

12.1 Replacing the O-rings

WARNING

Do not use lubricants such as grease to oil the O-rings.

Remove the worn O-rings and replace them with new ones. Carefully oil the new O-rings on the motor with MK-dent Premium Service Oil LU1011.

12.2 Replacing the light source


WARNING

Risk of burns due to hot lamp.

The light source generates a great deal of heat during operation.

→ Do not touch/replace the bulb until it has cooled down completely.


12.2.1 Replacing the light source for AM1116

Press the button (1) and pull off the sleeve (2) at the same time.

Turn the LED locking ring (3) until the LED lamp is free. Use a screwdriver to slide the LED lamp (4) out of the holder in the direction of the arrow. Place the new LED in the holder. Ensure that the contacts are in the correct position. Turn the LED locking ring (3) back again and refit the sleeve (2). If the lamp does not light up, turn the LED by 180°.

12.2.2 Replacing the light source for AM1006

Unscrew the LED cover in the direction of the arrow. Use a screwdriver to slide the LED lamp (1) out of the holder in the direction of the arrow. Place the new LED in the holder. Ensure that the contacts are in the correct position. Screw the LED cover (2) back on. If the lamp does not light up, turn the LED by 180°.

13 Technical data

Article no.	Drive speed [min ⁻¹]	Drive air pressure [bar]	Air consumption [l/min]	Return air pressure [bar]	Spray water pressure [bar]	Spray air pressure [bar]	Connection	Coating
AM0002	5.000 - 20.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	2-hole	Chrome
AM0004	5.000 - 20.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	4-hole	Chrome
AM0012	5.000 - 20.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	2-hole	Chrome
AM0014	5.000 - 20.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	4-hole	Chrome
AM1002	5.000 - 20.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	2-hole	Chrome
AM1004	5.000 - 20.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	4-hole	Chrome
AM1006	5.000 - 20.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	6-hole	Chrome
AM1116	5.000 - 22.000	2.2 – 2.8	max. 60	< 1,1	0,8 - 2,5	1,0 - 4,0	6-hole	Titanium

14 Explanation of symbols

	Attention – see section 1 Explanation of the warnings
	No sterilization
	No automatic processing / thermal disinfection
	Consult instructions for use and refer to the source of the electronic instructions for use
	CE marking with identification number of the notified body – indicates that the product is in compliance with European legislation for medical devices
	Medical device (CE conformity assessment procedure)
	Catalog number
	Serial number
	Manufacturer
	Date of manufacture
	Unique product identification – UDI (Unique Device Identifier) – indicates a barcode as containing Unique Device Identification
	HIBC Code (Health Industry Bar Code)

15 Guarantee and warranty

MK-dent assures you that your instrument has undergone a thorough **quality inspection** before leaving the factory. These quality inspections are carried out according to considerably stricter standards than the requirements of the product standard (ISO 14457 Dentistry - Handpieces and motors). In addition, we offer you comprehensive warranty services to ensure your satisfaction.

15.1 Guarantee

MK-dent grants the end customer a guarantee of perfect function and freedom from defects in material and workmanship from the date of purchase for a period of 24 months.

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapter, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (Instruments must be packed individually and with sufficient padding)
- No copy of the invoice for the device purchase is included (the invoice recipient must be a dental practice)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.2 Legal warranty

MK-dent grants the statutory warranty on all products to the contractual partner. This covers material and processing defects that were already present at the time of the transfer of risk or that occur within the statutory warranty period. Warranty claims will only be recognized if the product is submitted with proof of purchase in the form of an invoice or a copy of the delivery note (contractual partner, date of purchase, model and serial number must be clearly visible).

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (Instruments must be packed individually and with sufficient padding)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.3 Conditions

- Proof of purchase (invoice or delivery note) is required for guarantee or warranty claims.
- The guarantee or warranty is void if the product is tampered with or modified by unauthorized third parties.
- Repairs or replacements under warranty do not entitle the warranty period to be extended or restarted.
- MK-dent reserves the right to repair or replace defective instruments.

15.4 Utilization

In the event of a warranty claim, please contact MK-dent customer service directly. In the event of a warranty claim, please contact your contractual partner.

Contact MK-dent customer service

Phone: +49 4532 40049-0

Email: service@mk-dent.com

Please have the following information ready:

- Product model and serial number
- Proof of purchase with date
- Detailed description of the problem

Our customer service team will guide you through the warranty process and ensure a quick and satisfactory solution. After the repair, your instrument will be returned to you securely packaged and sealed together with a detailed test report.

15.5 Exclusion of liability

Within the scope of the guarantee or warranty, MK-dent is not liable for financial losses, downtimes, loan or rental equipment, loss of profit or similar consequences. The liability of MK-dent is limited to the purchase value of the product.

MK-dent is not liable for any defects and their consequences that may arise due to natural wear and tear, improper handling, improper cleaning or maintenance, non-compliance with operating or handling instructions, calcination or corrosion, contaminated air and water supply or chemical or electrical influences that are unusual or not permitted according to the instructions for use of MK-dent or other manufacturers' instructions for use.

Any liability is excluded if defects or their consequences are caused by interventions or modifications to the product by the end user or by a third party not authorized by MK-dent.

MK-dent reserves the right to change the warranty conditions at any time. Changes have no retroactive validity for products already purchased.

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