Instructions for use

MK-dent Propy Line Scaler

Version: 2025-03-17-RevA

Language Index

English Name	Native Name	ISO 639-1	Page
English	English	en	3
French	Français	fr	27
Spanish	Español	es	51
German	Deutsch	de	75
Italian	Italiano	it	100
Chinese	中文	zh	124
Swedish	Svenska	sv	149
Hungarian	Magyar	hu	173
Slovak	Slovenčina	sk	197
Czech	Česky	cs	221
Turkish	Türkçe	tr	245
Norwegian	Norsk	no	269
Bulgarian	Български	bg	293
Lithuanian	Lietuvių	lt	317
Portuguese (Portugal)	Português	pt	341
Arabic	العربية	ar	365
Greek	Αγγλικά	el	390

Table of contents

1	About these instructions for use	5
1.1	Explanation of warnings	5
1.2	Compliance	5
2	Intended purpose.....	6
2.1	Intended use.....	6
2.2	Indications	6
2.3	Contraindications	6
2.4	Intended patient groups.....	6
2.5	Intended users.....	6
2.6	User environment.....	6
3	Safety and requirements	7
3.1	General safety regulations	7
3.2	Reporting obligations	7
3.3	Side effects	8
3.4	Risk patient groups.....	8
4	Product description	8
4.1	Description	8
4.2	Combination with other products	8
4.3	Product components	9
4.4	Scope of delivery	9
4.5	Required materials that are not included	9
4.6	Configuration.....	10
5	Accessories, consumables and spare parts.....	11
5.1	Accessories	11
5.2	Consumables	11
5.3	Spare parts	11
6	Installation.....	11
6.1	Requirements for the environment	11
6.2	Unpacking.....	11
6.3	Assembly	11
6.4	Preparation for use.....	11
6.4.1	Inserting the scaler tip	12
6.4.2	Inserting the endo tip	13
6.4.3	Attaching the instrument.....	13
6.5	Startup and operation	13

7	Treatment	14
7.1	Patient, operator and assistant protection	14
7.2	Removing the scaler tip after use	14
7.3	Removing the endo tip after use	15
7.4	Scaler handpiece disassembly	15
8	Processing according to ISO 17664-1.....	16
8.1	Pre-treatment at the point of use	16
8.2	Manual processing	17
8.2.1	Manual cleaning and disinfection	18
8.2.2	Manual drying	18
8.3	Automatic processing	19
8.3.1	Automatic cleaning	19
8.3.2	Automatic disinfection and rinsing	19
8.3.3	Drying.....	19
8.4	Sterilization in an autoclave	20
8.4.1	Packaging	20
8.4.2	Sterilization in a steam sterilizer according to EN 13060 / EN ISO 17665	20
8.4.3	Sterilization parameters.....	21
8.4.4	After sterilization	21
8.5	Storage	21
9	Maintenance.....	22
9.1	Regular maintenance and checks	22
9.2	Cleaning in case of problems.....	22
10	Transport and disposal	22
11	Troubleshooting.....	22
12	Repair	23
13	Technical data.....	23
14	Explanation of symbols	24
15	Guarantee and warranty	25
15.1	Guarantee	25
15.2	Legal warranty	25
15.3	Conditions.....	25
15.4	Utilization.....	26
15.5	Exclusion of liability	26

1 About these instructions for use

These instructions for use apply to all MK-dent Prophy Line Scaler (see section **4 Product description**).

All products listed here are only intended for dental treatment in the field of dentistry. Any kind of misuse or modification of the product is not permitted and may lead to hazards. All products listed here may only be used by qualified users in accordance with **2.5 Intended users** for the intended purpose and in compliance with the instructions for use.

1.1 Explanation of warnings

These instructions for use are intended to describe the safe and effective operation of the product to the user. Before operating the product, you must read these instructions for use, observe all safety instructions and warnings and strictly adhere to them. Pay particular attention to all information and procedures described in section **3 Safety and requirements**.

These instructions for use are part of the product and must be kept within reach of the product so that they are always accessible.

To avoid personal injury and damage to property, observe the warnings and safety information contained in this document. These are specially marked:

WARNING

A WARNING alerts you to a potentially serious outcome, unwanted result or safety risk. Failure to observe a warning can result in the death or serious injury of the user or patient.

CAUTION

CAUTION alerts you to situations requiring particular care for the safe and effective use of the product. Disregard for precautionary measures can result in minor or moderately severe personal injuries or damage to the product or other property and may potentially lead to a remote risk of serious injuries and/or environmental pollution.

NOTE

NOTE

This identifies special tips, such as help to the user or for the improvement of a work process.

1.2 Compliance

This medical device fulfils the requirements of Directive 93/42/EEC.

If you have any further questions regarding the applicable national or international standards, please contact the manufacturer:

MK-dent GmbH, Marie-Curie-Str. 2, 22941 Bargteheide,
Germany
+49 4532 40049-0
info@mk-dent.com

2 Intended purpose

2.1 Intended use

The Prophy Line Scaler is an ultrasonic handpiece for the “drive” of scaler tips that are used for the removal of sub- and supragingival calculus during professional tooth cleaning by healthcare professionals. Ultrasonic handpieces and scaler tips are also used for cavity preparation, for cleaning, rinsing and preparing root canals, for loosening inlays, crowns and bridges and for periodontal therapy.

2.2 Indications

The Prophy Line Scaler is an ultrasonic handpiece for “driving” scaler tips that are used to remove sub- and supragingival calculus during professional tooth cleaning. Scaler tips are also used for cavity preparation, for cleaning, rinsing and preparing root canals, for loosening inlays, crowns and bridges and for periodontal therapy.

2.3 Contraindications

Due to the potential for electromagnetic interference, ultrasonic scalers are never used on patients with pacemakers, internal defibrillators, intracorporeal fluid pumps or other implantable electronic devices or in close proximity to sensitive patient monitoring devices such as pulse oximeters.

2.4 Intended patient groups

The target group is all patients who require sub- and/or supragingival scaling as part of prophylactic treatment. However, the decision to use sub- and/or supragingival scaling is the responsibility of the dentist or dental hygienist.

2.5 Intended users

The Prophy Line Scaler is intended for dentists or dental hygienists with appropriate training.

2.6 User environment

Application takes place in dental practices or dental facilities. Users wear protective gloves. The peripherals used meet the requirements of ISO 7494.

3 Safety and requirements

3.1 General safety regulations

WARNING

- This product is supplied non-sterile. Prepare it before first use in accordance with section **8 Processing according to ISO 17664-1**.
- For proper reprocessing, reprocess the device after each patient in accordance with section **8 Processing according to ISO 17664-1**. Safe and effective reprocessing after multiple uses without processing the device leads to mixed infections and prevents proper reprocessing due to atypical accumulation of contaminants, occlusions and cleaning agent residues – especially in cavities and supply lines.
- Do not use the ultrasonic scaler to process ceramic or metal dentures, as the high frequencies can cause damage.
- The scaler handpiece must not be used in conjunction with gaseous anesthetics in operating rooms.
- Functional check before each use. If you notice damage, irregular running noises, strong vibrations, unusual heating or other defects → notify Service
- Always use sufficient water during treatment to prevent the tip of the appliance from heating up. Otherwise, serious damage to the appliance, tooth, tooth root or gums may occur. To prevent the appliance from heating up, a water flow rate of at least 20 ml/min is required. Further details can be found in the instructions for use of the individual scaler attachments.
- This product must not be used in the vicinity of anyone wearing a pacemaker or suffering from hemophilia. Caution is also advised in the vicinity of pregnant women and children.
- This product may only be used by qualified personnel and only for dental treatment.
- Never disassemble or modify this product.
- This instrument must not be repaired by the user in the event of a malfunction. Contact your dealer or the MK-dent repair service.
- Please do not use the product if the performance of the product deteriorates or in case of malfunctions.
- Do not use this product if it has not been serviced properly.
- Observe the precautionary measures for use in order to ensure correct and safe functioning.

CAUTION

- Observe the instructions for the manufacturer's unit (e.g., Prophyl unit, dental chair, stand-alone unit): power setting for sub- and supragingival treatment.

3.2 Reporting obligations

In the event of serious deterioration in the patient's health during the use of this product, please report this to your specialized dealer or to us as the manufacturer – and to the national competent authority.

Our contact details are included in these instructions for use. The contact details of the notified body will also be provided to you upon request.

3.3 Side effects

There are currently no known side effects.

3.4 Risk patient groups

Never use this product on patients with:

- Pacemakers
- Internal defibrillators
- Intracorporeal fluid pumps
- Other implantable electronic devices
- Anyone suffering from endocarditis

4 Product description

4.1 Description

The MK-dent Prophy Line Scaler ultrasonic handpiece uses the piezo principle to generate 32,000 linear movements per second. The precisely defined movements of the attached scaler tip (must be purchased separately) oscillate back and forth between the tooth and gum in a controlled manner, thus protecting the adjacent tissue and the underlying nerves from injury. An LED ring, consisting of 6 LEDs, enables color-fast, glare-free and uniform illumination of the preparation site and the entire oral cavity (not with all models).

The MK-dent Prophy Line Scaler ultrasonic handpieces are compatible with the connections of manufacturers EMS, Acteon-Satelec and KaVo.

EMS, Acteon-Satelec and KaVo are registered trademarks. MK-dent has no economic connection with the aforementioned companies.

4.2 Combination with other products

CAUTION


Only use scaler tips that are compatible with the instrument to prevent damage to the instrument.

The MK-dent Prophy Line Scaler ultrasonic handpieces are compatible with the following scaler tips.

REF MK-dent Scaler	EMS	Acteon-Satelec
SC21E	X	---
SC21EL	X	---
SC21E2	X	---
SC21S	---	X
SC21SL	---	X
SC21KL	X	---

EMS and Acteon-Satelec are registered trademarks. MK-dent has no economic connection with the aforementioned companies.

4.3 Product components

4.4 Scope of delivery

The scope of delivery of the following components and accessories is included in the delivery unit:

REF MK-dent Scaler	OR3051 (1 pc.)	OR3052 (2 pcs.)	OR3053 (1 pc.)	OR3054 (1 pc.)	OR3055 (2 pcs.)
SC21E	X	X	---	---	---
SC21EL	X	---	---	---	---
SC21E2	X	---	---	---	---
SC21S	X	---	---	---	X
SC21SL	X	X	X	---	---
SC21KL	---	---	---	X	---

When you open the delivery package for the first time, make sure that the seal on the packaging is undamaged and that all components are included in the scope of delivery. If anything is missing, contact your dealer immediately.

4.5 Required materials that are not included

NOTE

NOTE

No scaler tip is included in the scope of delivery.
The scope of delivery does not include a dental treatment unit.

4.6 Configuration

The Prophy Line Scaler ultrasonic handpieces are connected to a dental treatment unit. There are different standardized connections depending on the manufacturer of the dental unit. The Prophy Line Scaler ultrasonic handpieces are compatible with the following tube connections:

REF	Variant EMS	Connection
SC21EL / SC21E2	Prophy Line ultrasonic scaler with LED	 EN060A / EN-061
SC21E	Prophy Line ultrasonic scaler	 EN-041

REF	Variant KaVo	Connection
SC21KL	Prophy Line ultrasonic scaler with LED	 REF 1.007.3995

REF	Variant Acteon-Satelec	Connection
SC21SL	Prophy Line ultrasonic scaler with LED	 F 12 609
SC21S	Prophy Line ultrasonic scaler	 F 12 281

EMS, Acteon-Satelec and KaVo are registered trademarks. MK-dent has no economic connection with the aforementioned companies.

5 Accessories, consumables and spare parts

5.1 Accessories

Not applicable

5.2 Consumables

Not applicable

5.3 Spare parts

Not applicable

6 Installation

6.1 Requirements for the environment

For regular maintenance of your dental treatment unit, please refer to the manufacturer's instructions for use. Ensure that the dental treatment unit is maintained accordingly.

6.2 Unpacking

When opening the package for the first time, make sure that the content according to **4.4 Scope of delivery** is enclosed with the product. Make sure that the Propy Line Scaler does not have any defects or damage such as cracks.

6.3 Assembly

Not applicable

6.4 Preparation for use

WARNING

Risk of infection for patients and users due to non-sterile instruments

→ Before initial use and after each treatment, the instrument and accessories must be processed in accordance with these instructions; see section **8**

Processing according to ISO 17664-1.

Malfunction and/or danger of infection due to moist and contaminated compressed air and contaminated water for cooling

→ The product may only be used in the field of dentistry with a dental treatment unit that has a CE marking or is otherwise certified/approved according to the valid national legal regulations and that is provided with a suitable supply of water and compressed air in compliance with the applicable standards.

CAUTION

A damaged instrument or damaged components can injure the patient, the user and third parties.

→ Test for proper functioning before each use and inspect the instrument visually for external damage.

→ Do not use the instrument and components if there are any abnormalities.

6.4.1 Inserting the scaler tip

WARNING

The tooth may be damaged if the scaler tip hits the tooth vertically during treatment.


CAUTION

The foot control of the treatment center must not be switched on while the scaler tip is being inserted.

Only use tips that are compatible with the scaler handpiece; see **4.2 Combination with other products**. When inserting the scaler tip, only use the tool from the corresponding scaler tip manufacturer to avoid overtightening the thread.

Use the appropriate tip template provided by the manufacturer to ensure that the scaler tip is in good enough condition before it is used. If not, replace it with a new one.

Screw the desired scaler tip into the handpiece before attaching it to the tube.

- Place the scaler tip on the thread at the top of the scaler handpiece (1).
- Turn the scaler tip clockwise and tighten it firmly (2).
- Remove the assembly tool. The instrument is now ready for use.
- Hold the handpiece between your thumb and index finger to work accurately. The treatment center can now be switched on.
- Set the vibration intensity and the water supply according to **13 Technical data**. For further information, refer to the corresponding instructions for use for the scaler tips.


6.4.2 Inserting the endo tip

WARNING


Do not operate the handpiece until the endo tip has been inserted into the root canal.

Only set the vibrations at a low intensity, as the tooth may otherwise be damaged. Do not apply too much pressure.

- Screw the file holder hand-tight onto the thread.
- Tighten the file holder using the appropriate manufacturer's tool.
- Insert the file into the file holder and tighten it with the tool.
- Set the vibration intensity and the water supply according to **13 Technical data**. For further information, please refer to the corresponding instructions for use for the endo tips.
- Now insert the file into the root canal and start the treatment.

6.4.3 Attaching the instrument

- Screw the desired scaler tip or endo tip into the handpiece before attaching it to the tube.
- Place the instrument exactly on the hose connection. Ensure that the contacts are in the correct position.
- Check that the instrument is securely attached to the hose by pulling gently.
- To remove, hold the instrument firmly and pull it off the tube sleeve.

6.5 Startup and operation

NOTE

NOTE

Before the start of each working day, the water supply systems must be flushed for at least 2 minutes without any instruments in place. Also rinse them for 20 to 30 seconds after each patient.

For startup and operation, follow the instructions for use provided by the manufacturer of your treatment unit and the scaler tip manufacturer.

7 Treatment

7.1 Patient, operator and assistant protection

WARNING

Before the first contact with the patient, the product must be properly processed to prevent infection. See section **8 Processing according to ISO 17664-1**.

Protect the patient accordingly before use.

During the treatment, the operator and assistants must wear gloves, safety goggles and masks to prevent infections.

CAUTION

Always stop working whenever operating noises become louder or irregular, there are stronger vibrations or in any case of overheating or visual damage. In these cases, contact your dealer or the MK-dent repair service.

7.2 Removing the scaler tip after use

CAUTION


The foot control of the treatment center must not be switched on while the scaler tip is being removed.

NOTE

NOTE


Only use the appropriate manufacturer's tool to remove the scaler tip in order to avoid overtightening the thread.

- Rinse the scaler for at least 30 seconds after each scaling process to remove coarse contaminants from the tip and handpiece.


- Place the tool on the scaler tip without touching it (1).
- Turn the scaler tip counterclockwise to loosen it from the thread (2).
- Pull off the scaler tip with the tool (3).

7.3 Removing the endo tip after use

- To remove the file, place the tool on the file holder and turn it counterclockwise to loosen the connection.
- Then remove the file.
- To remove the endo tip, place the tool on the thread and turn it counterclockwise to loosen the connection.
- Then remove the endo tip.

7.4 Scaler handpiece disassembly

- Unscrew the cap counterclockwise from the handpiece.
- Pull to remove the LED illuminated ring (only for scaler handpieces with a light) and the blue sealing ring.

8 Processing according to ISO 17664-1

NOTE

NOTE

The possible number of applications does not depend on the processing cycles but rather on the use of the device and the resulting wear and tear. The product's service life is therefore determined by signs of wear and tear and damage occurring during use.

Dispose of damaged products in accordance with section **10 Transport and disposal**. Do not use any products with obvious signs of wear or damage.

Use only cleaning agents / disinfectants that are authorized for the processing of medical devices. Follow the manufacturer's instructions for the cleaning agents / disinfectants.

If unsuitable cleaning agents or disinfectants are used, this can have negative effects on the instruments:

- Damage or corrosion
- Discoloration of the product
- Corrosion of metal parts
- Reduced service life

Follow the instructions of the cleaning agent / disinfectant manufacturer. Do not use any abrasive cleaning agents / disinfectants.

The manufacturer recommends only using cleaning agents / disinfectants compliant with EN ISO 15883-1 and -2 for mechanical cleaning / disinfecting.

The instruments can be sterilized. Do not use any sterilization temperatures above 135°C (275 °F).

Ultrasound is not recommended for processing the instruments, as this could reduce the product's service life.

8.1 Pre-treatment at the point of use

WARNING

Risk of infection from non-sterile instruments. Always wear protective gloves.

CAUTION

The use of disinfectant baths and disinfectants containing chlorine can cause defects and malfunctions of the instrument.

The cleaning result is diminished by surface drying of residues.

→ Flush the instrument immediately when you finish working on the patient.

→ Immediately remove any cement, blood or composite residues.

Manual pre-cleaning is also recommended in the case of automatic processing, such as to remove blockages in the spray channels.

- Pre-cleaning of the instruments should be carried out immediately after treatment, taking personal protection into account. The aim is to prevent organic material and chemical residues from drying in the cavity or on the outer parts of the instruments and to avoid contaminating the environment.
- Rinse the inner channels of the product still connected to the treatment center for at least 30 seconds.
- Remove the scaler tip or endodontic tip and remove the instrument from the treatment center.
- Remove adhering contaminants, body fluids and tissue with a disposable cloth / paper towel.
- Rinse the product with water after use (temperature below 35 °C (95 °F)).
- Ensure that the instrument is stored safely and transported in dry condition in a closed container to the processing site.
- Clean the instrument within one hour after each treatment to prevent it from drying out.
- Only use non-fixing cleaning agents or water with a temperature below 35 °C (95 °F), as there may otherwise be increased adhesion of contaminants, making further processing steps more difficult.

8.2 Manual processing

CAUTION

The use of manual pre-cleaning requires training of the employees who carry out this processing instruction. This is necessary to ensure successful cleaning.

Insufficient cleaning of the surface and water lines

→ Thorough manual pre-cleaning must be carried out to ensure successful cleaning.

Malfunctions can occur if cleaning is done in the ultrasonic device

→ Only clean the instrument as described below.

NOTE

NOTE

Use alcohol-based instrument disinfectant with a short residence time that does not cause any relevant adhesion of proteins on the handpieces. The instrument disinfectant must be bactericidal, fungicidal and antiviral according to RKI requirements. Never use disinfectants containing chlorine, as these can damage the handpieces. Do not use aldehydic disinfectants, as they are incompatible with the material. Only use disinfectants whose microbiological efficacy is confirmed or guaranteed by the manufacturer (e.g., VAH/DGHM registration or CE marking).

According to the RKI recommendation, thermal disinfection in a steam sterilizer is required following the manual cleaning and disinfecting of medical devices of the "semi-critical B" class, especially if the instrument disinfectant used for the manual disinfection does not simultaneously have bactericidal, fungicidal and antiviral efficacy. Alternatively, the manual cleaning can serve as preparation for mechanical cleaning.

Tools:

- Drinking water¹; min. 100 ml, 30 °C ± 5 °C (77 °F – 95 °F)
- Brush or soft toothbrush
- 1 % Dr. Weigert ersher® MediClean forte
- 70% isopropanol
- Disposable syringe

¹ Drinking water quality according to the EU Drinking Water Directive (total germ count max. 100 CfU/ml)

8.2.1 Manual cleaning and disinfection

- Remove the scaler tip or the endodontic tip from the handpiece using the respective assembly tool (see **7.2 Removing the scaler tip after use** or **7.3 Removing the endodontic tip after use**) and remove the instrument from the treatment unit.
- Brush the instrument under running drinking water¹ using a disposable toothbrush until no more contaminants are visible. Rinse the spray hole and the thread thoroughly.
- Then place the instrument in a cleaning solution (dosage 1%, Dr Weigert Neodisher® MediClean forte) and leave to soak for 10 minutes.
- Then rinse the instrument with demineralized water to remove the cleaning solution.
- Disassemble the instrument as described in **7.4 Scaler handpiece disassembly**.
- Place the disassembled handpiece in the alcohol-based disinfectant solution provided (70% isopropanol).
- Flush cavities with the disinfectant using a syringe.
- Leave the instruments to soak in the solution for 10 minutes.
- Then rinse all parts with demineralized water to remove the disinfectant solution.

8.2.2 Manual drying

Blow off the inside and outside of the instrument with medical compressed air until no more water droplets are visible.

¹ Drinking water quality according to the EU Drinking Water Directive (total germ count max. 100 CfU/ml)

8.3 Automatic processing

8.3.1 Automatic cleaning

NOTE

NOTE

MK-dent recommends thermal disinfection devices that satisfy the requirements of EN ISO 15883-1 and are approved by the manufacturer for dental handpieces. Cleaning should be performed at min. 55 °C (131 °F) for at least 6 minutes and disinfection at min. 90 °C (194 °F) for at least 5 minutes (for A0 value > 3000). A mildly alkaline cleaning agent with pH 9 to 11, such as Dr. Weigert Neodisher® MediClean forte is recommended.

CAUTION

For automatic cleaning of the interior, use the adapters specified by the device manufacturer.

Only process and disinfect instruments that are dry and free of residue using an automatic cleaning and disinfection process.

The Automatic cleaning process was validated using the WD BHT INNOVA® M3 device in program no. 14 (Cleaning) with Dr. Weigert Neodisher® MediClean forte 0.5 % as the cleaning agent. The instruments were connected by a silicone hose with a Luer-Lock adapter to the machine.

Parameters for automatic cleaning:

- Clean for 1 minute with 0.5% Neodisher® MediClean forte detergent at 30 °C (86 °F)
- Clean for 6 minutes with 0.5% Neodisher® MediClean forte at 55 °C (131 °F)
- Rinse for 1 minute with demineralized water

8.3.2 Automatic disinfection and rinsing

The automatic disinfection process was validated using the WD BHT INNOVA® M3 device in program no. 03 (Thermal disinfection) with Dr. Weigert Neodisher® MediClean forte 0.5 % as the cleaning agent. The instruments were connected by a silicone hose with a Luer-Lock adapter to the machine.

Parameters for automatic disinfection:

- Clean for 1 minute with 0.5% Neodisher® MediClean forte at 30 °C (86 °F)
- Clean for 6 minutes with 0.5% Neodisher® MediClean forte at 55 °C (131 °F)
- Rinse for 1 minute with demineralized water
- Rinse for 1 minute with demineralized water
- Carry out thermal disinfection for 5 minutes at 90 °C (194 °F) using low-germ (< 10 CFU/ml) and low-endotoxin (< 0.25 endotoxin units/ml), demineralized water (for an A0 value > 3,000).

8.3.3 Drying

The drying process is normally part of the cleaning process of your thermal disinfectant. Please follow the instructions for your thermal disinfectant. Allow the instrument to cool down before removing it from the adapter to ensure that this is safe and easier to do. If there is still visible contamination in the thermal disinfectant after treatment, the process must be repeated. To avoid any kind of damage to the instrument, make sure that the instrument is dry inside and outside and free of residue after each cycle. If moisture remains on or in the instrument after the cleaning process, dry it with a lint-free cloth. Also make sure that the cavities of the instruments are dry and dry them with medical compressed air if necessary.

8.4 Sterilization in an autoclave

8.4.1 Packaging

NOTE

NOTE

Use a validated packaging procedure that is suitable for sterilization. The sterilization bag must be big enough for the instrument so that the packaging does not get stretched. The sterilization packaging must conform to EN ISO 11607-1 and be suitable for the sterilization process. Pack each instrument individually in a simple sterilization pouch. Seal the sterilization pouch carefully and check the seal for leaks.

8.4.2 Sterilization in a steam sterilizer according to EN 13060 / EN ISO 17665

WARNING

Sterilization must be performed before the first use and after each patient to prevent direct and cross-infection.

To avoid infection, wear suitable protective clothing such as a mask, gloves and safety goggles.

Contact corrosion from moisture can cause damage to the product.

→ Take the instrument out of the steam sterilizer directly after the sterilization cycle.

This medical device can withstand temperatures up to max. 135 °C (275 °F).

CAUTION

Depending on the conditions of use of your sterilizer or your sterilization methods, this product may undergo discoloration.

NOTE

NOTE

The sterilizer manufacturer's instructions for operation and load configuration must be followed explicitly. When sterilizing multiple instruments in one sterilization cycle, make sure that the instruments do not touch each other to ensure the best possible steam penetration. Observe the information in the sterilizer manufacturer's instructions for use.

Put the instrument sealed in the sterilization bag on a sterilization tray to avoid direct contact with the inner chamber surface of the sterilizer and then place it into the autoclave. MK-dent recommends the use of a sterilizer in accordance with EN 13060, EN ISO 17665 or EN 285.

8.4.3 Sterilization parameters

Only use the following sterilization parameters to achieve an SAL of 10^{-6} :

Autoclave with a fractionated pre-vacuum (recommended): at least 4 minutes at min. 132 °C (270 °F).

The sterilization was validated with the W&H LISA 517 device in a fractionated pre-vacuum process at 132 °C (270 °F) for 4 minutes. According to the sterilizer manufacturer's instructions for use, the drying time for the sterilization process is calculated by the device's software.

The instructions listed above for processing a medical device for reuse have been validated and found to be suitable by the medical device manufacturer. The user is responsible for ensuring that the processing actually performed, including the devices and materials used in the processing facility and any personnel involved, will achieve the intended result. This requires verification and/or validation and routine monitoring of the process.

8.4.4 After sterilization

When the sterilization cycle is finished, take the instrument out of the autoclave. When doing so, be sure to use appropriate personal protective equipment, as the instruments will be hot.

8.5 Storage

After sterilization and drying, check that there are no drops of residual water. Store the processed instruments in a dry, dark and cool room, protected from dust and with the lowest possible microbial load.

Temperature	+5 °C - +40 °C (+41 °F - +104 °F)
Relative humidity	< 80 %, non-condensing
Air pressure	500 hPa - 1060 hPa (7.2 psi - 15 psi)

CAUTION

Instrument malfunction during operation after heavily cooled storage conditions
→ Allow heavily cooled instruments to warm up to room temperature before starting to use them.

The instrument must be completely dried before storage. Any remaining water can result in discoloration or rust.

Do not store the instrument in a place where it could be exposed to water.

Do not store the instrument in places where chemicals are stored or where corrosive gases exist.

9 Maintenance

9.1 Regular maintenance and checks

WARNING

Inappropriate maintenance or service can result in malfunctions and/or deterioration.

→ Always conduct proper maintenance and servicing.

CAUTION

For safe use of the product, perform checks and maintenance of the instrument according to these instructions for use.

9.2 Cleaning in case of problems

In case of problems, clean each part of the instrument. The instrument must be disconnected from the hose for cleaning. Instructions can be found in section **8.1 Pre-treatment at the point of use**.

10 Transport and disposal

CAUTION

Users or third parties may be infected by contaminated medical devices. Process the product before transport or disposal or dispose of it properly in compliance with the applicable regulations for the disposal of contaminated medical waste.

Process the instrument in accordance with section **8 Processing in accordance with ISO 17664-1** and return it to us in the original packaging. We will take care of the professional disposal of the device and the original packaging for you.

11 Troubleshooting

NOTE

NOTE

This instrument must not be repaired by the user in the event of a malfunction. Please contact your dealer or the MK-dent repair service.

12 Repair

CAUTION

Contaminated medical devices can infect users or third parties. Process the instrument before dispatching it.

Process this instrument according to section **8 Processing according to ISO 17664-1** and send it to your dealer or directly to MK-dent if it is defective. Please contact MK-dent customer service to arrange for shipping or collection for factory repair. Questions about noise, maintenance, processing and functionality can often be clarified in a conversation, making it unnecessary to send in the product.

13 Technical data

	SC21E, SC21EL, SC21E2 (EMS connection)	SC21S, SC21SL (Acteon-Satelec connection)	SC21KL (KaVo connection)
Electronic voltage input [V]	190 - 230	190 - 230	190 - 230
Range of movement of the instrument [°]	1 ~ 100	1 ~ 100	1 ~ 100
Movement frequency of the instrument [kHz]	28 ± 3	28 ± 3	28 ± 3
Output power [W]	3 - 20	3 - 20	3 - 20
Spray water pressure [bar]	0.1 - 5.0	0.1 - 5.0	0.1 - 5.0
Spray water consumption [ml/min] according to EN ISO 18397	> 20	> 20	> 20

EMS, Acteon-Satelec and KaVo are registered trademarks. MK-dent has no economic connection with the aforementioned companies.

14 Explanation of symbols

	Attention – see section 1 Explanation of the warnings
	Sterilizable in a steam sterilizer (autoclave) up to the specified temperature
	Suitable for automatic processing in thermal disinfection devices
	Consult instructions for use and refer to the source of the electronic instructions for use
	CE marking with identification number of the notified body – indicates that the product is in compliance with European legislation for medical devices
	Medical device (CE conformity assessment procedure)
	Catalog number
	Serial number
	Legal manufacturer
	Date of manufacture
	Unique product identification – UDI (Unique Device Identifier) – indicates a barcode as containing Unique Device Identification
	HIBC Code (Health Industry Bar Code)

15 Guarantee and warranty

MK-dent assures you that your instrument has undergone a thorough **quality inspection** before leaving the factory. These quality inspections are carried out according to considerably stricter standards than the requirements of the product standard (ISO 18397 Dentistry – Powered scaler). In addition, we offer you comprehensive warranty services to ensure your satisfaction.

15.1 Guarantee

MK-dent grants the end customer a guarantee of perfect function and freedom from defects in material and workmanship from the date of purchase for a period of 12 months.

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (Instruments must be packed individually and with sufficient padding)
- No enclosure of the invoice copy of the device purchase (invoice recipient must be a dental practice)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.2 Legal warranty

MK-dent grants the statutory warranty on all products to the contractual partner. This covers material and processing defects that were already present at the time of the transfer of risk or that occur within the statutory warranty period. Warranty claims will only be recognized if the product is submitted with proof of purchase in the form of an invoice or a copy of the delivery note (contractual partner, date of purchase, model and serial number must be clearly visible).

Exclusions:

- Wear parts, e.g. collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (instruments must be packed individually and with sufficient padding)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.3 Conditions

- Proof of purchase (invoice or delivery note) is required for guarantee or warranty claims.
- The guarantee or warranty is void if the product is tampered with or modified by unauthorized third parties.
- Repairs or replacements under warranty do not entitle the warranty period to be extended or restarted.
- MK-dent reserves the right to repair or replace defective instruments.

15.4 Utilization

In the event of a warranty claim, please contact MK-dent customer service directly. In the event of a warranty claim, please contact your contractual partner.

Contact MK-dent customer service

Phone: +49 4532 40049-0

Email: service@mk-dent.com

Please have the following information ready:

- Product model and serial number
- Proof of purchase with date
- Detailed description of the problem

Our customer service team will guide you through the warranty process and ensure a quick and satisfactory solution. After the repair, your instrument will be returned to you securely packaged and sealed together with a detailed test report.

15.5 Exclusion of liability

Within the scope of the guarantee or warranty, MK-dent is not liable for financial losses, downtimes, loan or rental equipment, loss of profit or similar consequences. The liability of MK-dent is limited to the purchase value of the product.

MK-dent is not liable for any defects and their consequences that may arise due to natural wear and tear, improper handling, improper cleaning or maintenance, non-compliance with operating or handling instructions, calcination or corrosion, contaminated air and water supply or chemical or electrical influences that are unusual or not permitted according to the instructions for use of MK-dent or other manufacturers' instructions for use.

Any liability is excluded if defects or their consequences are caused by interventions or modifications to the product by the end user or by a third party not authorized by MK-dent.

MK-dent reserves the right to change the warranty conditions at any time. Changes have no retroactive validity for products already purchased.

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