Instructions for use

MK-dent Quick Couplings

Version: 2025-03-19-RevA

Language Index

English Name	Native Name	ISO 639-1	Page
English	English	en	3
French	Français	fr	25
Spanish	Español	es	47
German	Deutsch	de	69
Italian	Italiano	it	91
Chinese	中文	zh	113
Swedish	Svenska	sv	136
Hungarian	Magyar	hu	158
Slovak	Slovenčina	sk	180
Czech	Česky	cs	202
Turkish	Türkçe	tr	224
Norwegian	Norsk	no	246
Bulgarian	Български	bg	277
Lithuanian	Lietuvių	lt	290
Portuguese (Portugal)	Português	pt	312
Arabic	العربية	ar	334
Greek	Αγγλικά	el	356

Table of contents

1	About these instructions for use	5
1.1	Explanation of warnings	5
1.2	Compliance	5
2	Intended purpose.....	6
2.1	Intended use.....	6
2.2	Indications	6
2.3	Contraindications	6
2.4	Intended patient groups.....	6
2.5	Intended users.....	6
2.6	User environment.....	6
3	Safety and requirements	7
3.1	General safety regulations	7
3.2	Reporting obligations	7
3.3	Side effects	7
3.4	Risk patient groups.....	7
4	Product description	8
4.1	Description	8
4.2	Combination with other products	8
4.3	Product components	8
4.4	Scope of delivery	8
4.5	Required materials that are not included	9
4.6	Configurations	9
5	Accessories, consumables and spare parts.....	9
5.1	Accessories	9
5.2	Consumables	9
5.3	Spare parts	10
6	Installation.....	11
6.1	Requirements for the environment	11
6.2	Unpacking.....	11
6.3	Assembly	11
6.4	Preparing the connection	11
6.4.1	Connection to the handpiece.....	11
6.4.2	Connection to the hose coupling	11
6.5	Startup and operation	12
6.5.1	Check of correct operation	12

6.5.2	Setting the compressed air	12
6.5.3	Water regulation	12
7	Treatment	13
7.1	Patient, operator and assistant protection	13
7.2	Disconnecting after use	13
8	Processing according to ISO 17664-1	14
8.1	Pre-treatment at the point of use	14
8.2	Manual processing	15
8.2.1	Manual cleaning of the surface	15
8.2.2	Manual interior cleaning and drying	15
8.2.3	Manual disinfection	16
8.3	Automatic processing	16
8.4	Sterilization	16
8.5	Storage	17
9	Maintenance	17
9.1	Regular maintenance and checks	17
10	Transport and disposal	18
11	Troubleshooting	18
12	Repair	18
12.1	Replacing the O-rings	18
12.2	Replacing the light source	19
12.2.1	Turbines with hose connection	19
12.2.2	Turbines with coupling connection	19
13	Technical data	20
14	Explanation of symbols	22
15	Guarantee and warranty	23
15.1	Guarantee	23
15.2	Legal warranty	23
15.3	Conditions	23
15.4	Utilization	24
15.5	Exclusion of liability	24

1 About these instructions for use

These instructions for use apply to all MK-dent quick couplings (see section **4 Product description**).

All products listed here are only intended for dental treatment in the field of dentistry. Any kind of misuse or modification of the product is not permitted and may lead to hazards. All products listed here may only be used by qualified users in accordance with **2.5 Intended users** for the intended purpose and in compliance with the instructions for use.

1.1 Explanation of warnings

These instructions for use are intended to describe the safe and effective operation of the product to the user. Before operating the product, you must read these instructions for use, observe all safety instructions and warnings and strictly adhere to them. Pay particular attention to all information and procedures described in section **3 Safety and requirements**.

These instructions for use are part of the product and must be kept within reach of the product so that they are always accessible.

To avoid personal injury and damage to property, observe the warnings and safety information contained in this document. These are specially marked:

WARNING

A WARNING alerts you to a potentially serious outcome, unwanted result or safety risk. Failure to observe a warning can result in the death or serious injury of the user or patient.

CAUTION

CAUTION alerts you to situations requiring particular care for the safe and effective use of the product. Disregard for precautionary measures can result in minor or moderately severe personal injuries or damage to the product or other property and may potentially lead to a remote risk of serious injuries and/or environmental pollution.

NOTE

NOTE

This identifies special tips, such as help to the user or for the improvement of a work process.

1.2 Compliance

This medical device fulfills the requirements of Directive 93/42/EEC.

If you have any further questions regarding the applicable national or international standards, please contact the manufacturer:

MK-dent GmbH, Marie-Curie-Str. 2, 22941 Bargteheide,
Germany
+49 4532 40049-0
info@mk-dent.com

2 Intended purpose

2.1 Intended use

Couplings are an accessory for turbines and are designed to connect the hose from the dental unit to the turbine to conduct air, water and light. Couplings are medical devices without direct patient contact.

2.2 Indications

Not applicable. Couplings are an accessory for turbines without direct patient contact.

2.3 Contraindications

There are currently no known contraindications.

2.4 Intended patient groups

Not applicable. Couplings are an accessory for dental turbines without direct patient contact.

2.5 Intended users

A coupling is only intended for use by licensed dentists or dental hygienists with appropriate training. Secondary users are the practice staff who are responsible for transport, processing and cleaning.

2.6 User environment

Application takes place in dental practices or dental facilities. Users wear protective gloves. The peripherals used fulfill the requirements of ISO 7494.

3 Safety and requirements

3.1 General safety regulations

WARNING

- This product is supplied non-sterile. Disinfect it before first use in accordance with section **8 Processing according to ISO 17664-1**.
- For proper reprocessing, reprocess the product after each patient in accordance with section **8 Reprocessing according to ISO 17664-1**. Safe and effective reprocessing after multiple uses without reprocessing the device leads to mixed infections and prevents proper reprocessing due to atypical accumulation of contaminants, occlusions and cleaning agent residues – especially in cavities and supply lines.
- This product may only be used by qualified personnel and only for air-operated handpieces.
- Never disassemble or modify this product.
- This instrument must not be repaired by the user in the event of a malfunction. Contact your dealer or the MK-dent repair service.
- Please do not use the product if the performance of the product deteriorates or if there are any malfunctions.
- Do not use this product if it has not been serviced properly.
- Observe the precautionary measures for use in order to ensure correct and safe functioning.

CAUTION

- Before initial use and after each treatment, the instrument and accessories must be prepared in accordance with the instructions for use.
- The product may only be used in the field of dentistry with a dental treatment unit that has a CE marking or is otherwise certified/approved according to the valid national legal regulations and that is provided with a suitable supply of water and compressed air in compliance with the applicable standards.
- Before each use, carry out a function test and check the instrument visually for external damage and contamination.

3.2 Reporting obligations

In the event of a serious deterioration in the patient's health during the use of this product, please report this to your specialized dealer or to us as the manufacturer – and also to the national competent authority.

Our contact details are included in these operating instructions. The contact details of the notified body will also be provided to you upon request.

3.3 Side effects

There are currently no known side effects.

3.4 Risk patient groups

There are currently no known risk patient groups.

4 Product description

4.1 Description

Couplings are an accessory for air-operated handpieces and are designed to connect the hose from the dental unit to the handpiece in order to conduct air, water and, if necessary, light.

Thanks to the quick-coupling system, replacing air-operated handpieces is easy, enabling immediate use. No complicated installation is required, allowing for simple integration in your treatment system.

Driven by compressed air and water from the treatment unit, a coupling neither influences other medical devices nor is it affected by electromagnetic waves.

4.2 Combination with other products

REF Coupling	QC4012K, QC4014K, QC5015K, QC5016K, QC5016KW, QC5116K, QC5116KW, QC6016K, QC6016KW, QC6116K, QC6116KW	QC4014NT, QC5016NT, QC6016NT	QC5016SW, QC6016SW	QC4014W, QC5016W
REF MK- dent Product	HE20K/KL, HE21K/KL, HE22K/KL, HB21K/KL, HB23K/KL, HP21KL, HC20K/KL, HC21K/KL, HC22K/KL, HP22KL, PR1011K	HE20N/NL, HE21N/NL, HE22N/NL, HP22NL, HP21NL, HE21NL, HE22NL, HE21N, HE22N, PR1011N	HE21SL, HP21SL, PR1011S	HE21W/WL, HE22W/WL, HP21WL, HP22WL, PR1011W
Connection	KaVo	NSK	Sirona	W&H

KaVo, Sirona, W&H and NSK are registered trademarks. MK-dent has no economic connection with the aforementioned companies.

4.3 Product components

Not applicable

4.4 Scope of delivery

The scope of delivery of the following components and accessories is included in the delivery unit:

Product group	Key	O-rings
Couplings with KaVo connection	KE1029	3 x OR4011BK, 2 x OR4011WK
Couplings with NSK connection	---	---
Couplings with Sirona connection	KE1029S	4 x OR4011S
Couplings with W&H connection	---	2 x OR4011BW, 1 x OR4012BW

KaVo, Sirona, W&H and NSK are registered trademarks. MK-dent has no economic connection with the aforementioned companies.

When you open the delivery package for the first time, make sure that the seal on the packaging is undamaged and that all components are included in the scope of delivery. If anything is missing, contact your dealer immediately.

4.5 Required materials that are not included

NOTE

NOTE

The scope of delivery does not include an air-operated handpiece.

The scope of delivery does not include MK-dent Premium Service Oil LU1011.

4.6 Configurations

Not applicable

5 Accessories, consumables and spare parts

5.1 Accessories

Not applicable

5.2 Consumables

Not applicable

5.3 Spare parts

CAUTION

The use of unauthorized accessories or unauthorized modifications to the product can lead to injuries or malfunctions of the instrument.

→ Only use spare parts and components that have been approved by MK-dent for combination with the product.

QC4014NT, QC5016NT, QC6016NT	
Replacement coupling with LED/XENON/without light	QC6016NT / QC5016NT / QC4014NT
Replacement LED / XENON	BU8012N / BU7012N

QC6016SW, QC5016SW	
Replacement coupling with LED / XENON	QC6016SW / QC5016SW
Replacement LED / XENON	BU8012SB1 / BU7012SB
Replacement O-ring	OR4011S
O-ring set 10x OR4011S	OR4011S/10

QC6016KW, QC6116KW / QC5016KW, QC5116KW	
Replacement coupling with LED/XENON	QC6016KW, QC6116KW / QC5016KW, QC5116KW
Replacement LED / XENON	BU8012 / BU7012
Replacement O-ring white	OR4011WK
Replacement O-ring black	OR4011BK
O-ring set 3x OR4011BK, 2x OR4011WK	OR4011BWK

QC4014W, QC5016W	
Replacement coupling without / with power connection	QC4014W / QC5016W
Replacement O-ring small	OR4012BW
Replacement O-ring large	OR4011BW
O-ring set 2x OR4011BW, 1x OR4012BW	OR4011BBW

QC6016K, QC6116K / QC5016K, QC5116K / QC4012K, QC4014K, QC5015K	
Replacement coupling with LED / XENON	QC6016K, QC6116K / QC5016K, QC5116K
Replacement coupling with light guide	QC5015K
Replacement coupling without light	QC4012K, QC4014K
Replacement LED / XENON	BU8012 / BU7012
Replacement O-ring white	OR4011WK
Replacement O-ring black	OR4011BK
O-ring set 3x OR4011BK, 2x OR4011WK	OR4011BWK

6 Installation

6.1 Requirements for the environment

To supply the coupling with clean and dry air, clean the compressor air filter in your dental practice. Do not supply moist compressed air. For drainage and regular maintenance of your compressor, please refer to its instructions for use. Ensure that the compressor is maintained accordingly.

6.2 Unpacking

When opening the package for the first time, make sure that the content according to **4.4 Scope of delivery** is enclosed with the product. Make sure that the coupling does not have any defects or damage such as cracks.

6.3 Assembly

Not applicable

6.4 Preparing the connection

WARNING

When connecting your air-operated handpiece to this product, check that the connection is secure. An unstable connection during a treatment can lead to injuries due to the body or other parts falling.

Never pull the handpiece out of your coupling while the handpiece is in operation. This can lead to injuries.

Check that the product is free of damage, cracks, breakage and corrosion and that it is safe to use.

CAUTION

The couplings must be disinfected before initial use and after each treatment. Only use dry, uncontaminated compressed air. Dirty or moist compressed air leads to premature wear.

When an air turbine, handpiece, etc. is removed, water and oil remain at the coupling connection. Remove water and oil completely. Make sure that no water or oil remains in the coupling of your turbine hose. If any is present, remove it with a dry cloth.

6.4.1 Connection to the handpiece

To connect to the handpiece, please refer to the instrument manufacturer's instructions for use.

6.4.2 Connection to the hose coupling

Attach the coupling to the hose and tighten the union nut of the hose firmly.

6.5 Startup and operation

6.5.1 Check of correct operation

WARNING

Before starting any treatment, make sure that the coupling and instrument are correctly and firmly connected. Only use instruments that are safe to use.

Press the foot switch of your treatment unit once or twice for about one second and check – when air and water are sprayed – whether the water supply is sufficient and no leakage of air or water can be detected at the connection of the coupling.

If you have problems, read section **11 Troubleshooting**. Use the water volume regulator on your treatment unit to set the water supply.


6.5.2 Setting the compressed air

WARNING

The compressed air at the coupling must be set in accordance with the handpiece's instructions for use.

6.5.3 Water regulation

Turn in the direction of the arrow to regulate the water.

The maximum amount of spray water is reached when the markings on the rings are on top of each other.

7 Treatment

7.1 Patient, operator and assistant protection

WARNING

Before the first contact with the patient, the product must be properly processed to prevent infection. See section **8 Processing according to ISO 17664-1**.

Protect the patient accordingly before use.

During treatment, the operator and assistants must wear gloves, protective goggles and masks to prevent infection.

CAUTION

Always stop working whenever operating noises become louder or irregular, there are stronger vibrations or in any case of overheating or visual damage. In these cases, contact your dealer or the MK-dent repair service.

7.2 Disconnecting after use

WARNING

The product must not be dropped.

Disconnect the coupling in accordance with your handpiece's instructions for use. When reconnecting, first clean the coupling, turbine hose and the inside of the connection in accordance with **6.4 Preparing the connection**.

8 Processing according to ISO 17664-1

NOTE

NOTE

The possible number of applications does not depend on the processing cycles but rather on the use of the device and the resulting wear and tear. The product's service life is therefore determined by signs of wear and tear and damage occurring during use.

Dispose of damaged products in accordance with section **10 Transport and disposal**. Do not use any products with obvious signs of wear or damage.

Only use cleaning agents / disinfectants that are authorized for the processing of medical devices. Follow the manufacturer's instructions for the cleaning agents / disinfectants.

If unsuitable cleaning agents or disinfectants are used, this can have negative effects on the instruments:

- Damage or corrosion
- Discoloration of the product
- Corrosion of metal parts
- Reduced service life

Follow the instructions of the cleaning agent / disinfectant manufacturer. Only use non-fixing cleaning agents / disinfectants.

Couplings cannot be sterilized. Sterilization is not absolutely necessary for the intended use.

Ultrasonic cleaning is not recommended for processing the instruments, as this could reduce the product's service life.

8.1 Pre-treatment at the point of use

WARNING

Risk of infection from non-sterile instruments. Always wear protective gloves.

CAUTION

The use of disinfectant baths and disinfectants containing chlorine can cause defects and malfunctions of the instrument.

The cleaning result is diminished by surface drying of residues.

→ Flush the instrument immediately when you finish working on the patient.

Do not place the instrument in solutions or similar.

NOTE

NOTE

Clean the connection hose and the inside of the instrument connector according to **6.4 Preparation for use** before reconnecting the instrument.

- Pre-cleaning of the instruments should be carried out immediately after treatment, taking personal protection into account. The aim is to prevent organic material and chemical residues from drying on the coupling and to avoid contaminating the environment.
- Remove adhering contaminants with a disposable cloth / paper towel.
- Rinse the product with water after use (temperature below 35 °C (95 °F)).
- Only use non-fixing cleaning agents or water at a temperature below 35 °C (95 °F), as there may otherwise be increased adhesion of contaminants, making further processing steps more difficult.
- Clean the coupling within one hour after each treatment to prevent it from drying out.
- Ensure that the instrument is stored safely and transported in dry condition in a closed container to the processing site.

8.2 Manual processing

CAUTION

Insufficient cleaning of the surface and of water and air supply lines.

Thorough manual pre-cleaning must be carried out to ensure successful cleaning.

The use of manual pre-cleaning requires training of the employees who carry out this processing instruction. This is necessary to ensure successful cleaning.

Malfunctions can occur if cleaning is done in the ultrasonic machine.

→ Clean the instrument exclusively as described below.

8.2.1 Manual cleaning of the surface

Tools:

- Drinking water¹ 30 °C ± 5 °C (77 °F - 95 °F)
- Brush or medium-hard toothbrush

Brush the coupling under running drinking water¹.

8.2.2 Manual interior cleaning and drying

Flush the spray air and spray water line with drinking water¹.

Blow the coupling with medical compressed air from the outside and from the inside until no more drops of water are visible.

¹ Drinking water quality according to the EU Drinking Water Directive (total germ count max. 100 CfU/ml)
MK-dent-Instructions-for-Use-Quick-Couplings-Multilingual-RevA
Version: 2025-03-19-RevA

8.2.3 Manual disinfection

WARNING

Couplings cannot be sterilized.

The use of disinfection baths and disinfectants containing chlorine can lead to defects and malfunctions in the coupling.

→ Only clean and disinfect the coupling manually

NOTE

NOTE

Use alcohol-based instrument disinfectant with a short residence time that does not cause any relevant adhesion of proteins on the handpieces. The instrument disinfectant must be bactericidal, fungicidal and antiviral according to RKI requirements. Never use disinfectants containing chlorine, as these can damage the handpieces. Do not use aldehydic disinfectants, as they are incompatible with the material. Only use disinfectants whose microbiological efficacy is confirmed or guaranteed by the manufacturer (e.g., VAH/DGHM registration or CE marking).

Tools:

- Instrument disinfectant
- Disposable wipe

Spray the disinfectant onto a disposable cloth and then wipe the coupling. Do not place the instrument in disinfectant baths!

Observe the exposure time specified by the disinfectant manufacturer. Follow the disinfectant's instructions for use.

8.3 Automatic processing

Not applicable

WARNING

Automatic processing is not to be used for couplings.

8.4 Sterilization

Not applicable

WARNING

Sterilization in an autoclave is not to be used for couplings.

8.5 Storage

After disinfection and drying, check that there are no drops of residual water. Store the processed instruments in a dry, dark and cool room, protected from dust and with the lowest possible microbial load.

Temperature	-20 °C - +70 °C (-4 °F - +158 °F)
Relative humidity	5% - 95%, non-condensing
Air pressure	700 hPa - 1060 hPa (10 psi - 15 psi)

CAUTION

Instrument malfunction during operation after heavily cooled storage conditions
→ Allow heavily cooled instruments to warm up to room temperature before starting to use them.

The instrument must be completely dry before storage. Any remaining water can result in discoloration or rust.

Do not store the coupling in a place where it may be exposed to water.

Do not store the coupling in places where chemicals are stored or where corrosive gases exist.

9 Maintenance

9.1 Regular maintenance and checks

WARNING

Inappropriate maintenance or service can lead to malfunctions and/or deterioration.

→ Always conduct proper maintenance and servicing.

CAUTION

For safe use of the product, perform checks and maintenance of the instrument according to these instructions for use.

10 Transport and disposal

CAUTION

Users or third parties may be infected by contaminated medical devices. Process the product before transport or disposal or dispose of it properly in compliance with the applicable regulations for the disposal of contaminated medical waste.

Process the instrument in accordance with section **8 Processing in accordance with ISO 17664-1** and return it to us in its original packaging. We will take care of the professional disposal of the device and the original packaging for you.

11 Troubleshooting

CAUTION

This instrument must not be repaired by the user in the event of a malfunction. Please contact your dealer or the MK-dent repair service.

Process this instrument according to section **8 Processing according to ISO 17664-1** and send it to your dealer or directly to MK-dent if it is defective. Please contact MK-dent customer service to arrange for shipping or collection for factory repair. Questions about noise, maintenance, processing and functionality can often be clarified in a conversation, making it unnecessary to send in the product.

12 Repair

12.1 Replacing the O-rings

WARNING

Do not use lubricants such as grease to oil the O-rings.

Remove the worn O-rings and replace them with new ones. Carefully oil the new O-rings on the coupling with MK-dent Premium Service Oil LU1011.

12.2 Replacing the light source

WARNING

Risk of burns due to hot lamp.

The light source generates a great deal of heat during operation.


→ Do not touch/replace the bulb until it has cooled down completely.

NOTE


NOTE

LEDs work with an exact pole position. If the LED does not work after installation, install it again rotated by 180°.


12.2.1 Turbines with hose connection

Unscrew the threaded ring in the direction of the arrow.


Remove the inner adapter.


Remove the defective bulb and insert the new one into the bulb socket.

12.2.2 Turbines with coupling connection

Remove the lamp cover in the direction of the arrow.



Remove the defective bulb.



Insert a new bulb and replace the protective lamp cover.


13 Technical data

		QC4014NT	QC5016NT	QC6016NT
	Optimum luminosity at 184 mA	-	3.2 - 3.5 V	3.2 - 3.5 V
	Light	-	XENON	LED
	Type	4-hole	6-hole	
	Suitable connection	NSK Mach-Lite		
	Suitable replacement bulb	-	BU7012N	BU8012N
	Coating	Chrome		

		QC5016SW	QC6016SW
	Optimum luminosity at 184 mA	3.2 - 3.5 V	
	Light	XENON	LED
	Type	6-hole	
	Suitable connection	Sirona	
	Suitable replacement bulb	BU7012SB	BU8012SB1
	Coating	Chrome	

		QC5116KW	QC5016KW	QC6116KW	QC6016KW
	Optimum luminosity at 184 mA	3.2 - 3.5 V			
	Light	XENON		LED	
	Type	6-hole			
	Suitable connection	MK-dent, KaVo Multiflex			
	Suitable replacement bulb	BU7012		BU8012	
	Water regulation	Yes			
	Coating	Titanium	Chrome	Titanium	Chrome

		QC4012K	QC4014K	QC5015K	QC5116K	QC5016K	QC6116K	QC6016K	
	Optimum luminosity at 184 mA	-	-	-	3.2 - 3.5 V				
	Light	-	-	Glass fiber rod	XENON		LED		
	Type	2-hole	4-hole	5-hole	6-hole				
	Suitable connection	MK-dent, KaVo Multiflex							
	Suitable replacement bulb	-	-	-	BU7012		BU8012		
	Coating	Chrome			Titanium	Chrome	Titanium	Chrome	

		QC4014W	QC5016W
	Power supply	-	Yes
	Type	4-hole	6-hole
	Suitable connection	W&H Roto Quick	
	Coating	Chrome	Chrome

KaVo, Sirona, W&H and NSK are registered trademarks. MK-dent has no economic connection with the aforementioned companies.

14 Explanation of symbols

	Caution – see section 1 Explanation of the warnings
	No sterilization
	No automatic processing / thermal disinfection
	Consult instructions for use and refer to the source of the electronic instructions for use
	CE marking with identification number of the notified body – indicates that the product is in compliance with European legislation for medical devices
	Medical device (CE conformity assessment procedure)
	Catalog number
	Serial number
	Manufacturer
	Date of manufacture
	Unique product identification – UDI (Unique Device Identifier) – indicates a barcode as containing Unique Device Identification
	HIBC Code (Health Industry Bar Code)

15 Guarantee and warranty

MK-dent assures you that your instrument has undergone a thorough **quality inspection** before leaving the factory. These quality inspections are carried out according to considerably stricter standards than the requirements of the product standard (ISO 14457 Dentistry – Handpieces and motors). In addition, we offer you comprehensive warranty services to ensure your satisfaction.

15.1 Guarantee

MK-dent grants the end customer a guarantee of perfect function and freedom from defects in material and workmanship from the date of purchase for a period of 24 months.

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (Instruments must be packed individually and with sufficient padding)
- No enclosure of the invoice copy of the device purchase (invoice recipient must be a dental practice)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.2 Legal warranty

MK-dent grants the statutory warranty on all products to the contractual partner. This covers material and processing defects that were already present at the time of the transfer of risk or that occur within the statutory warranty period. Warranty claims will only be recognized if the product is submitted with proof of purchase in the form of an invoice or a copy of the delivery note (contractual partner, date of purchase, model and serial number must be clearly visible).

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (instruments must be packed individually and with sufficient padding)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.3 Conditions

- Proof of purchase (invoice or delivery note) is required for guarantee or warranty claims.
- The guarantee or warranty is void if the product is tampered with or modified by unauthorized third parties.
- Repairs or replacements under warranty do not entitle the warranty period to be extended or restarted.
- MK-dent reserves the right to repair or replace defective instruments.

15.4 Utilization

In the event of a warranty claim, please contact MK-dent customer service directly. In the event of a warranty claim, please contact your contractual partner.

Contact MK-dent customer service

Phone: +49 4532 40049-0

E-mail: service@mk-dent.com

Please have the following information ready:

- Product model and serial number
- Proof of purchase with date
- Detailed description of the problem

Our customer service team will guide you through the warranty process and ensure a quick and satisfactory solution. After the repair, your instrument will be returned to you securely packaged and sealed together with a detailed test report.

15.5 Exclusion of liability

Within the scope of the guarantee or warranty, MK-dent is not liable for financial losses, downtimes, loan or rental equipment, loss of profit or similar consequences. The liability of MK-dent is limited to the purchase value of the product.

MK-dent is not liable for any defects and their consequences that may arise due to natural wear and tear, improper handling, improper cleaning or maintenance, non-compliance with operating or handling instructions, calcination or corrosion, contaminated air and water supply or chemical or electrical influences that are unusual or not permitted according to the instructions for use of MK-dent or other manufacturers' instructions for use.

Any liability is excluded if defects or their consequences are caused by interventions or modifications to the product by the end user or by a third party not authorized by MK-dent.

MK-dent reserves the right to change the warranty conditions at any time. Changes have no retroactive validity for products already purchased.

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