Instructions for use

MK-dent straight and contra angle Handpieces

Version: 2025-12-23-RevB

Language Index

English Name	Native Name	ISO 639-1	Page
English	English	en	3
French	Français	fr	37
Spanish	Español	es	71
German	Deutsch	de	105
Italian	Italiano	it	140
Chinese	中文	zh	174
Swedish	Svenska	sv	208
Hungarian	Magyar	hu	242
Slovak	Slovenčina	sk	276
Czech	Česky	cs	310
Turkish	Türkçe	tr	344
Norwegian	Norsk	no	378
Bulgarian	Български	bg	412
Lithuanian	Lietuvių	lt	446
Portuguese (Portugal)	Português	pt	480
Arabic	العربية	ar	514
Greek	Αγγλικά	el	547

Table of contents

1	About these instructions for use.....	6
1.1	Explanation of warnings	6
1.2	Compliance	6
2	Intended purpose	7
2.1	Intended use	7
2.2	Indications.....	7
2.3	Contraindications.....	7
2.4	Intended patient groups	7
2.5	Intended users	7
2.6	User environment	7
3	Safety and requirements	8
3.1	General safety regulations	8
3.2	Reporting obligation	8
3.3	Side effects.....	8
3.4	Risk patient groups	8
4	Product description.....	9
4.1	Description.....	9
4.2	Combination with other products	10
4.3	Product components	10
4.4	Scope of delivery.....	11
4.5	Required materials that are not included.....	11
4.6	Configurations.....	11
5	Accessories, consumables and spare parts.....	12
5.1	Accessories	12
5.2	Consumables.....	13
5.3	Spare parts	13
6	Installation.....	14
6.1	Requirements for the environment.....	14
6.2	Unpacking	14
6.3	Assembly.....	14
6.4	Preparation for use	14
6.4.1	Connecting to the motor	15
6.5	Startup and operation	16
6.5.1	Before use	16
6.5.2	Checking for correct operation	16
6.5.3	Testing the water quantity	16
6.5.4	Use of rotating instruments	17

6.5.4.1	Contra-angle handpieces	18
6.5.4.2	Straight handpieces	18
6.5.4.3	Insert prophy cups or polishing brushes	18
7	Treatment.....	19
7.1	Patient, operator and assistant protection.....	19
7.2	Removing rotating instruments	19
7.2.1	Contra-angle handpieces.....	19
7.2.2	Straight handpieces.....	20
7.2.3	Removing prophy cups or polishing brushes	20
8	Processing according to ISO 17664-1	21
8.1	Pre-treatment at the point of use.....	22
8.1.1	Preparation for cleaning.....	22
8.1.2	Manual pre-cleaning of the spray nozzle	23
8.1.3	Manual pre-cleaning of the instrument	23
8.1.4	Manual pre-cleaning of the glass rod of handpieces with light	24
8.2	Manual processing.....	25
8.2.1	Manual cleaning	25
8.2.2	Manual disinfection, rinsing and drying.....	25
8.3	Automatic processing	26
8.3.1	Automatic cleaning, disinfection and rinsing	26
8.3.2	Automatic drying.....	26
8.3.3	Maintenance, checks and testing.....	27
8.4	Sterilization in an autoclave.....	27
8.4.1	Packaging.....	27
8.4.2	Sterilization in a steam sterilizer according to EN 13060 / EN ISO 17665	28
8.4.3	Sterilization parameters	28
8.4.4	After sterilization.....	29
8.5	Storage.....	29
9	Maintenance.....	29
9.1	Regular maintenance and checks	29
9.2	Cleaning in case of problems	29
10	Transport and disposal	30
11	Troubleshooting.....	30
12	Repair	31
13	Technical data.....	32
14	Explanation of symbols	34
15	Guarantee and warranty	35
15.1	Guarantee	35

15.2 Legal warranty.....35

15.3 Conditions 36

15.4 Utilization 36

15.5 Exclusion of liability.....36

1 About these instructions for use

These instructions for use apply to all MK-dent straight handpieces and contra-angle handpieces (see section **4 Product description**).

All products listed here are only intended for dental treatment in the field of dentistry. Any kind of misuse or modification of the product is not permitted and may lead to hazards. All products listed here may only be used by qualified users in accordance with **2.5 Intended users** for the intended purpose and in compliance with the instructions for use.

1.1 Explanation of warnings

These instructions for use are intended to describe the safe and effective operation of the product to the user. Before operating the product, you must read these instructions for use, observe all safety instructions and warnings and strictly adhere to them. Pay particular attention to all information and procedures described in section **3 Safety and requirements**.

These instructions for use are part of the product and must be kept within reach of the product so that they are always accessible.

To avoid personal injury and damage to property, observe the warnings and safety information contained in this document. These are specially marked:

WARNING

A WARNING alerts you to a potentially serious outcome, unwanted result or safety risk. Failure to observe a warning can result in the death or serious injury of the user or patient.

CAUTION

CAUTION alerts you to situations requiring particular care for the safe and effective use of the product. Disregard for precautionary measure can result in minor or moderately severe personal injuries or damages to the product or other property and may potentially lead to a remote risk of serious injuries and/or environmental pollution.

NOTE

NOTE

This identifies special tips, such as help to the user or for the improvement of a work process.

1.2 Compliance

If you have any further questions regarding the applicable national or international standards, please contact the manufacturer:

MK-dent GmbH, Marie-Curie-Str. 2, 22941 Bargteheide,
Germany
+49 4532 40049-0
info@mk-dent.com

2 Intended purpose

2.1 Intended use

Straight handpieces (slowly rotating, max. 40,000 min⁻¹): drives for rotating instruments for use in the preparation of dentures and for polishing tooth and restoration surfaces outside of the oral cavity.

Contra-angle handpieces (slowly rotating, 2,703 – 40,000 min⁻¹): drives for rotating instruments for polishing teeth and dentures and for preparing cavities.

Contra-angle handpieces (up to 200,000 min⁻¹): drives for rotating instruments for the removal of carious material, the preparation of cavities and crowns and the removal of fillings.

2.2 Indications

Diseases of the periodontium, teeth and jaws (ICD-10 codes).

Handpieces and contra angle handpieces are motor-driven preparation instruments with different rotation speed ranges. They use mechanical energy to rotate the drill head to change structures, such as tooth substances, jawbones or filling materials, whereby it is the responsibility of the user to decide which tissue is to be removed.

Depending on the work process, drills (e.g., rose-head drill, fissure drill), diamonds (actually: steel shafts fitted with diamond grains of different shapes and sizes), grinding bodies (corundum) or rubber polishers can be clamped into the handpiece.

For excavating (removal of carious dentine), a rotation speed between 1,000 min⁻¹ and 1,500 min⁻¹ is recommended; where necessary, the maximum rotation speed is selected for the preparation (grinding) of enamel and dentine for fitting a crown.

For dental technology work (corrections) on denture parts, a straight handpiece is preferable because it is more robust and the instruments clamped in have a longer HP (handpiece) shaft.

2.3 Contraindications

There are currently no known contraindications.

2.4 Intended patient groups

The decision to use a contra-angle handpiece, however, is the sole responsibility of the dentist. Depending on the assessment of the user, the products can be used on all patients who require dental treatment.

2.5 Intended users

MK-dent handpieces and contra-angle handpieces are intended exclusively for use by accredited dentists. Secondary users are the practice staff who are responsible for transport, processing and cleaning.

2.6 User environment

Application takes place in dental practices or dental facilities. Users wear protective gloves. The peripherals used fulfil the requirements of ISO 7494.

3 Safety and requirements

3.1 General safety regulations

WARNING

- This product is supplied non-sterile. Prepare it before first use in accordance with section **8 Processing according to ISO 17664-1**.
- For proper reprocessing, process the device after each patient in accordance with section **8 Reprocessing according to ISO 17664-1**. Safe and effective processing after multiple uses without processing the device leads to mixed infections and prevents proper processing due to atypical accumulation of contaminants, occlusions and cleaning agent residues – especially in cavities and supply lines.
- This product may be used only by qualified personnel and only for dental treatment.
- Never disassemble or modify this product.
- This instrument must not be repaired by the user in the event of a malfunction. Contact your dealer or the MK-dent repair service.
- Please do not use the product if the performance of the product deteriorates or in case of malfunctions.
- Do not use this product if it has not been serviced properly.
- Observe the precautionary measures for use in order to ensure correct and safe functioning.

CAUTION

- The product may only be used in the field of dentistry with a dental treatment unit that has a CE marking or is otherwise certified/approved according to the valid national legal regulations and that is provided with an appropriate water and compressed air supply in compliance with the applicable standards.
- Test for proper functioning (collet holding force, spray water flow rate and heat development) before each use and inspect the instrument visually for external damage and contamination.
- If the amount of spray water is too low, the instrument may overheat and cause injuries.

3.2 Reporting obligation

In the event of serious deterioration in the patient's health during the use of this product, please report this to your specialized dealer or to us as the manufacturer and to the national competent authority.

Our contact details are included in these instructions for use. The contact details of the notified body will also be provided to you upon request.

3.3 Side effects

There are currently no known side effects.

3.4 Risk patient groups

There are currently no known risk patient groups.

4 Product description

4.1 Description

MK-dent straight and contra-angle handpieces are available in the following variants: PRIME LINE, CLASSIC LINE, ECO LINE, BASIC LINE

Within these product lines, both straight (transmission 1:1 = 40,000 min⁻¹) and contra-angle handpieces are available (transmission 4:1 = 5,000 min⁻¹, 1:5 = 200,000 min⁻¹, plus additional variants generating different rotation speeds with the same input).

The compact shape of the handpieces makes switching between treatment units easy and requires very little space for storage.

No complicated installation is required, allowing for simple integration in your treatment system.

Article no.	Light	Straight handpiece	Contra-angle handpiece complete	Contra-angle handpiece body	Contra-angle handpiece head	Motor connection ISO 3964	HP drill	RA drill	FG drill	Prophy cup
LB01	--	X	--	--	--	X	X	--	--	--
LB01W	--	X	--	--	--	X	X	--	--	--
LB02	--	X	--	--	--	X	X	--	--	--
LB1111	--	--	X	--	--	X	--	X	--	--
LB1111C	--	--	X	--	--	X	--	X	--	--
LB11W	--	--	X	--	--	X	--	X	--	--
LB4111	--	--	X	--	--	X	--	X	--	--
LB4111P	--	--	X	--	--	X	--	--	--	X
LBH11	--	--	--	--	X	--	--	X	--	--
LBH11C	--	--	--	--	X	--	--	X	--	--
LBH11P	--	--	--	--	X	--	--	--	--	X
LE01	--	X	--	--	--	X	X	--	--	--
LE01L	X	X	--	--	--	X	X	--	--	--
LE11	--	--	X	--	--	X	--	X	--	--
LE11L	X	--	X	--	--	X	--	X	--	--
LE148	--	--	X	--	--	X	--	X	--	--
LE148L	X	--	X	--	--	X	--	X	--	--
LE148LP	X	--	X	--	--	X	--	--	--	X
LE148P	--	--	X	--	--	X	--	--	--	X
LE15	--	--	X	--	--	X	--	--	X	--
LE15L	X	--	X	--	--	X	--	--	X	--
LE16	--	--	X	--	--	X	--	--	X	--
LE16L	X	--	X	--	--	X	--	--	X	--
LE21	--	--	X	--	--	X	--	X	--	--
LE21L	X	--	X	--	--	X	--	X	--	--
LE21P	--	--	X	--	--	X	--	--	--	X
LE27	--	--	X	--	--	X	--	X	--	--
LE27L	X	--	X	--	--	X	--	X	--	--
LE54	--	--	X	--	--	X	--	X	--	--
LE54L	X	--	X	--	--	X	--	X	--	--
LE54LP	X	--	X	--	--	X	--	--	--	X
LE54P	--	--	X	--	--	X	--	--	--	X
LE74	--	--	X	--	--	X	--	X	--	--
LE74L	X	--	X	--	--	X	--	X	--	--
LC01	--	X	--	--	--	X	X	--	--	--
LC01L	X	X	--	--	--	X	X	--	--	--
LC11	--	--	--	X	--	X	--	--	--	--

Article no.	Light	Straight handpiece	Contra-angle handpiece complete	Contra-angle handpiece body	Contra-angle handpiece head	Motor connection ISO 3964	HP drill	RA drill	FG drill	Prophy cup
LC11L	X	--	--	X	--	X	--	--	--	--
LC27	--	--	--	X	--	X	--	--	--	--
LC27L	X	--	--	X	--	X	--	--	--	--
LC74	--	--	--	X	--	X	--	--	--	--
LC74L	X	--	--	X	--	X	--	--	--	--
LCH11	--	--	--	--	X	--	--	X	--	--
LCH16	--	--	--	--	X	--	--	--	X	--
LCH21	--	--	--	--	X	--	--	X	--	--
LCH21P	--	--	--	--	X	--	--	--	--	X
LP01L	X	X	--	--	--	X	X	--	--	--
LP11L	X	--	X	--	--	X	--	X	--	--
LP15L	X	--	X	--	--	X	--	--	X	--
LP27L	X	--	X	--	--	X	--	X	--	--
LP54L	X	--	X	--	--	X	--	X	--	--
LP74L	X	--	X	--	--	X	--	X	--	--
LPH11	--	--	--	--	X	--	--	X	--	--
LPH16	--	--	--	--	X	--	--	--	X	--
LPH21	--	--	--	--	X	--	--	X	--	--

4.2 Combination with other products

REF MK-dent straight or contra-angle handpiece	LB01, LB01W, LB02, LB1111, LB1111C, LB11W, LB4111, LB4111P, LE01, LE01L, LE11, LE11L, LE148, LE148L, LE148LP, LE148P, LE15, LE15L, LE16, LE16L, LE21, LE21L, LE21P, LE27, LE27L, LE54, LE54L, LE54LP, LE54P, LE74, LE74L, LC01, LC01L, LC11, LC11L, LC27, LC27L, LC74, LC74L, LP01L, LP11L, LP15L, LP27L, LP54L, LP74L
REF MK-dent air motor (20,000 min ⁻¹)	AM0002, AM0004, AM0012, AM0014, AM1002, AM1004, AM1006, AM1116
REF MK-dent electric motor (40,000 min ⁻¹)	EM1190K
Connection*	ISO connection*

*The straight and contra-angle handpieces are compatible with all motors with a connection according to DIN EN ISO 3964.

4.3 Product components

Not applicable

4.4 Scope of delivery

The scope of delivery of the following components and accessories is included in the delivery unit:

Product variant	AC0001 Nozzle pin	AC0003 Dental brush for maintenance of the spindle	LT1013P Oil adapter for the motor connection according to DIN EN ISO 3964	LT1014P Oil adapter for maintenance of the spindle
Contra-angle handpiece, complete	X	X	X	X
Contra-angle body / straight handpiece	X	X	X	--
Contra-angle head	--	X	X	X

When you open the delivery package for the first time, make sure that the seal on the packaging is undamaged and that all components are included in the scope of delivery. If anything is missing, contact your dealer immediately.

4.5 Required materials that are not included

NOTE

NOTE

The scope of delivery does not include a motor.

The scope of delivery does not include an instrument holder, item no. RT2050.

The scope of delivery does not include a holding force measuring tool, item no. RT1020.


The scope of delivery does not include MK-dent Premium Service Oil LU1011.

4.6 Configurations

The configurations of the straight handpieces and contra-angle handpieces are described in **4.1 Description** and **5.3 Spare parts**.

5 Accessories, consumables and spare parts

5.1 Accessories

Holding force measuring tool RT1020	
Instrument holder RT2050	
MK-dent Premium Service Oil LU1011	

5.2 Consumables

Not applicable

5.3 Spare parts

CAUTION

The use of unauthorized accessories or unauthorized modifications to the product can lead to injuries or instrument malfunctions.

→ Only use spare parts and components that have been approved by MK-dent for combination with the product.

	Straight handpiece
Prime Line replacement handpiece	LP01L
Classic Line replacement handpiece	LC01, LC01L
Eco Line replacement handpiece	LE01, LE01L
Basic Line replacement handpiece	LB01, LB01W, LB02

	Contra-angle handpiece, complete	Replacement body	Replacement head
Prime Line replacement contra-angle handpiece	LP11L LP15L LP27L LP54L LP74L	-- -- -- -- --	LPH11 -- LPH11 LPH21 LPH11
Classic Line replacement contra-angle handpiece	--	LC11, LC11L, LC27, LC27L LC74, LC74L	LCH11, LCH16, LCH21, LCH21P, LPH11, LPH16, LPH21
Eco Line replacement contra-angle handpiece	LE11/L LE15/L LE16/L LE21/L LE21P LE27/L LE54/L LE54P/LP LE74/L LE148/L LE148P/LP	-- -- -- -- -- -- -- -- -- -- --	LCH11 -- LCH16 LCH21 LCH21P LCH11 LCH21 LCH21P LCH11 LCH21 LCH21P
Basic Line replacement contra-angle handpiece	LB1111 LB1111C LB11W LB4111 LB4111P	-- -- -- -- --	LBH11 LBH11C -- LBH11 LBH11P

6 Installation

6.1 Requirements for the environment

To supply the products with clean and dry air, clean the compressor air filter in your dental practice. Do not supply moist compressed air. For drainage and regular maintenance of your compressor, please read the operating instructions for your air compressor. Ensure that the compressor is maintained in accordance with its instructions for use.

6.2 Unpacking

When opening the package for the first time, make sure that the content according to **4.4 Scope of delivery** is enclosed with the product. Make sure that the straight or contra-angle handpiece has no defects or damage such as cracks.

6.3 Assembly

Not applicable

6.4 Preparation for use

WARNING

Danger of an infection for patients and users due to unsterile instruments
 → Before initial use and after each treatment, the instrument and accessories must be processed according to these instruction for use; see section **8 Processing according to ISO 17664-1**.

Malfunction and/or danger of infection from moist and contaminated compressed air and contaminated water for cooling
 → The product may only be used in the field of dentistry with a dental treatment unit that has a CE marking or is otherwise certified/approved according to the valid national legal regulations and that is provided with a suitable supply of water and compressed air in compliance with the applicable standards.

CAUTION

A damaged instrument or damaged components can injure the patient, the user and third parties.
 → Test for proper functioning before each use and inspect the instrument visually for external damages.
 → Do not use the instrument and components if there are any abnormalities.

If a straight or contra-angle handpiece is removed, water and oil will remain on the motor connection. Remove the water and oil completely. Confirm that no water or oil remains in the coupling of your hose. If there is still any residue left, remove it with a dry cloth.

Make sure that no water or oil remains in the connection of the instrument. Remove residues with a cotton swab.

6.4.1 Connecting to the motor


WARNING

Malfunction due to decoupling of the instrument during treatment.
An instrument that is not properly engaged can come loose from the motor during treatment.
→ Check that the instrument is securely attached to the motor by pulling on it before each treatment.

CAUTION

Damage caused by pulling off and putting on the instrument while the motor is rotating.
→ Damage to the driving pin
→ Never put on or take off the body while the motor is rotating

- Place the instrument exactly onto the motor coupling and press it down firmly until you can hear the connection snap in place.
- On handpieces with light, set the correct position of the light conductor by turning the instrument firmly until the lock for the correct seat of the light guide snaps in.
- Check that the instrument is securely attached to the motor coupling by pulling on it.
- For removal, hold the instrument and pull it off by turning it carefully.

6.5 Startup and operation

6.5.1 Before use

NOTE

NOTE

Before the start of every working day, the systems conducting water must be flushed for at least 2 minutes without any instruments attached. Another flushing process must be performed for 20 to 30 seconds after each patient.

6.5.2 Checking for correct operation

WARNING

Check the safe connection to the motor after you have heard a locking sound. An instable connection leads to injuries caused by the instrument or other parts falling off.

Check for the absence of damage, cracks, breakage and any corrosion on the product and that it is safe to use.

If you have problems, read section **11 Troubleshooting**. Use the water volume regulator on your treatment center to set the water supply.

6.5.3 Testing the water quantity

WARNING

Injury to the pulp and tooth.

A low quantity of spray water can result in overheating of the instrument and cause injuries.

→ Set the water quantity for the spray cooling to min. 50 ml/min.

Check the water quantity with a measuring cup. Check the spray water tubes. If the water output is too low or a nozzle is clogged, clean the spray nozzle with the enclosed nozzle pin. Also see **8.1.2 Manual pre-cleaning of the spray nozzle**.

6.5.4 Use of rotating instruments

WARNING

risk of injury to patients if the instrument falls out.

→ After each use of rotating instruments, check the holding force of the collet by pulling on the rotating instrument. Please be sure to do this while the actuator is stopped.

Risk of injury to patients caused by damage to the instrument. Use of rotating instruments other than those listed below.

→ Use only rotating instruments with specifications that are equivalent to the data listed below.

→ Observe the instructions for use and proper use of the rotating instruments.

Risk of injury due to the use of worn rotating instruments. Rotating instruments can fall out during the treatment and injure the patient.

→ Do not use rotating instruments with worn shafts.

Risk of injury to users due to rotating instruments. Cuts and potential infections can be caused as a result.

→ Use gloves or finger protection


NOTE

NOTE

Use only rotating instruments made of steel or tungsten carbide that comply with DIN EN ISO 1797 and meet the following criteria:

	Shaft diameter			Maximum clamping length		Snap-on prophyl cups
	Type 1 (CA, RA): 2.334 to 2.350 mm with groove	Type 2 (HP): 2.334 to 2.350 mm cylindrical	Type 3 (FG): 1.590 to 1.600 mm cylindrical	12 mm	30 mm	
REF MK-dent straight or angle handpiece						
LB1111, LB1111C, LB111W, LB4111, LBH11, LBH11C, LE11, LE11L, LE148, LE148L, LE21, LE21L, LE27, LE27L, LE54, LE54L, LE74, LE74L, LCH11, LCH21, LP11L, LP27L, LP54L, LP74L, LPH11, LPH21	X	--	--	X	--	--
LB01, LB01W, LB02, LE01, LE01L, LC01, LC01L, LP01L	--	X	--	--	X	--
LE15, LE15L, LE16, LE16L, LCH16, LP15L, LPH16	--	--	X	X	--	--
LB4111P, LBH11P, LE148P, LE148LP, LE21P, LE54P, LE54LP, LCH21P	--	--	--	--	--	X

6.5.4.1 Contra-angle handpieces


- Press the button firmly with your thumb and simultaneously insert the rotating instrument with short rotating movements until it stops.
- Pull on the drill to check that it sits firmly.
- Check the holding force of the spindle at least once a week using the MK-dent holding force measuring tool RT1020 or another suitable device according to ISO 8325. The holding force of the spindle must not fall below 2.5 kg. MK-dent recommends replacing the spindle when the holding force is below 3 kg.

6.5.4.2 Straight handpieces

- Turn the tension ring in the direction of the arrow as far as it will go and insert the rotating instrument into the spindle.
- Turn the tension ring back to the starting position.
- Pull on the rotating instrument to check that it is secure.

6.5.4.3 Insert prophyl cups or polishing brushes

Snap-on and screw-type polishing bits (prophy cups and brushes) according to DIN EN ISO 13295 can be used.

- Screw the bit into the head or slide it onto the head.
- Check that the bit is securely in place.

7 Treatment

7.1 Patient, operator and assistant protection

WARNING

Before the first contact with the patient, the product must be properly processed to prevent infection. See section **8 Processing according to ISO 17664-1**.

Protect the patient appropriately before use.

During the treatment, the operator and assistants must wear gloves, safety goggles and masks to prevent infections.

CAUTION

Stop working whenever operating noises get louder or irregular, if there are stronger vibrations, or in any case of overheating or visual damage. In such cases, contact your dealer or the MK-dent repair service.

MK-dent recommends lubricating the instrument after a treatment duration of 30 minutes, using MK-dent Premium Service Oil LU1011; see section **8.3.3 Maintenance, checks and testing**. If the instrument is to be used again directly after it is lubricated, it must run in idle for 15 - 20 seconds to remove excess oil. If the instrument is to be used again later on the same patient, retain it in an instrument holder (item no. RT2050) so that any excess oil can drain.

7.2 Removing rotating instruments


WARNING

Risk of injury from rotating instruments. Cuts and damage to the clamping system.

- Never touch rotating instruments.
- Never press the button while the instrument is rotating.
- Remove the rotating instrument after each treatment to prevent injuries and infections when putting it down

7.2.1 Contra-angle handpieces


- Use your thumb to press the button firmly after the rotating instrument has come to a complete stop and pull it out at the same time.

7.2.2 Straight handpieces

- After the rotating instrument has come to a complete stop, turn the tension ring in the direction of the arrow up to the stop and take out the rotating instrument.
- Never lock the handpiece without rotating the instrument or test pin.

7.2.3 Removing prophyl cups or polishing brushes

- Screw the bit out of the head or take it off.

8 Processing according to ISO 17664-1

NOTE

NOTE

The lifecycle of the medical device is validated for 250 processing cycles. After 250 processing cycles, no significant negative effect on the product condition from processing was able to be detected. The possible number of applications does not depend on the processing cycles but rather on the application of the device and the wear and tear caused by this. The product's service life is therefore determined by signs of wear and damage occurring during use.

Dispose of damaged products in accordance with section **10 Transport and disposal**. Do not use any products with obvious signs of wear or damage.

Only use cleaning agents / disinfectants that are authorized for the processing of medical devices. Follow the manufacturer's instructions for the cleaning agents / disinfectants.

If unsuitable cleaning agents or disinfectants are used, this can have negative effects on the instruments:

- Damage or corrosion
- Discoloration of the product
- Corrosion of metal parts
- Reduced service life

Follow the instructions of the cleaning agent / disinfectant manufacturer. Do not use any abrasive cleaning agents / disinfectants.

The manufacturer recommends only using cleaning agents / disinfectants compliant with EN ISO 15883-1 and -2 for automated cleaning/disinfecting.

The instruments can be sterilized. Do not use any sterilization temperatures above 135 °C (275 °F).

Ultrasonic cleaning is not recommended for processing the instruments, as the ball bearings can be damaged, thereby reducing the product's service life.

8.1 Pre-treatment at the point of use

WARNING

Risk of infection from non-sterile instruments. Always wear protective gloves.

CAUTION

Using disinfecting baths and disinfectants containing chlorine can cause defects and malfunctions of the instrument.

The cleaning result is diminished by surface-drying of residues.

→ Flush the transfer instrument immediately when you finish working on the patient.

→ Immediately remove any cement, blood or composite residues.

Do not place the instrument in solutions or similar.

NOTE

NOTE

Clean the connection hose and on the inside of the instrument connector according to **6.4 Preparation for use** before reconnecting the instrument.

- Flush the inner tubes of the product according to the RKI specifications for 20 seconds while it is still connected to the treatment unit.
- The instruments should be pre-cleaned directly after the treatment for the sake of personal safety. The aim is to prevent organic material and chemical residues from drying in the cavity or on the outer parts of the instruments and to avoid contaminating the environment.
- Remove adhering contaminants, body fluids and tissue using a disposable cloth/paper tissue.
- Rinse the product with water after use (temperature below 35 °C (95 °F)).
- Do not use any abrasive cleaning agents or water with temperatures above 35 °C (95 °F), as this can lead to increasing adhesion of the contaminants, making the further processing steps more difficult.
- Clean the instrument within one hour after each treatment to prevent surface drying.
- Ensure that the instrument is stored safely and transported in dry condition in a closed container to the place of processing.

8.1.1 Preparation for cleaning

Remove the drill bit from the collet and pull the instrument off the motor or quick coupling. Take off all removable parts (if any) and wash each part under running water.

Manual pre-cleaning is also recommended in the case of automatic processing, such as for removing blockages in the spray channels.


8.1.2 Manual pre-cleaning of the spray nozzle

WARNING

Injury to the pulp and tooth. A low quantity of spray water can result in overheating of the instrument and corresponding injuries.

→ Check spray water tubes and, if necessary, clean the spray nozzles with the nozzle pin (AC0001) included in the delivery.

Carefully remove clogging in the spray tubes with the enclosed nozzle pin.

8.1.3 Manual pre-cleaning of the instrument

CAUTION

Insufficient cleaning of the surface and water and air tubes.

Thorough manual pre-cleaning must be carried out to ensure successful cleaning.

The use of manual pre-cleaning requires training of the employees who carry out this processing instruction. This is necessary to ensure successful cleaning.

Tools:

- Drinking water¹ min. 100 ml 30°C ± 5°C (77 °F – 95 °F)
- Brush or soft toothbrush
- Dental brush (AC0003)
- 50 ml syringe or cannula
- Soft lint-free cloth

CAUTION


Risk of infection from dental brushes reused multiple times

→ To avoid cross infection, use a new dental brush for each instrument that is to be cleaned

Malfunctions can occur if cleaning is done in the ultrasonic machine

→ Clean the instrument only as described below

¹ Drinking water quality according to the EU Drinking Water Directive (total germ count max. 100 CfU/ml)
MK-dent-Instructions-for-Use-Low-Speed-Handpieces-Multilingual-RevB
Version: 2025-12-23-RevB

- Carefully brush the complete instrument under running drinking water¹ using a brush or soft toothbrush to pre-clean the instrument.
- To pre-clean the spindle, insert the enclosed dental brush (AC0003) into the collet and carefully loosen/remove any deposits and contaminants in rotating movements.
- Repeat this process several times and flush the dental brush with drinking water¹ each time.
- To pre-clean the water and air tubes, apply the syringe to the bottom end of the handpiece and flush out the hollow spaces at least five times.

NOTE**NOTE**

Repeat these steps several times in case of heavy contamination.

8.1.4 Manual pre-cleaning of the glass rod of handpieces with light

**CAUTION**

Risk of scratching and defects on the glass rod of the instruments if not handled with care.

→ Be extremely careful when you clean the glass rod

Conduct a visual inspection after each cleaning of the glass rod. Do not use the handpiece or angle handpiece if the glass rod is damaged and contact your dealer or the MK-dent repair service.

Wash the glass rod with warm drinking water¹ using a soft lint-free glass cleaning cloth.

Dry the glass rod using a compressed air gun or do so carefully using a soft lint-free glass cleaning cloth.

¹ Drinking water quality according to the EU Drinking Water Directive (total germ count max. 100 CfU/ml)

8.2 Manual processing

NOTE

NOTE

Use alcohol-based instrument disinfectant with a short residence time that does not cause any relevant adhesion of proteins on the handpieces. The instrument disinfectant must be bactericidal, fungicidal and antiviral according to RKI requirements. Never use disinfectants containing chlorine, as these can damage the handpieces. Do not use aldehydic disinfectants, as they are incompatible with the material. Only use disinfectants whose microbiological efficacy is confirmed or guaranteed by the manufacturer (e.g., VAH/DGHM registration or CE marking).

According to the RKI recommendation, thermal disinfection in a steam sterilizer is required following the manual cleaning and disinfecting of medical devices of the "semi-critical B" class, especially if the instrument disinfectant used for the manual disinfection does not simultaneously have bactericidal, fungicidal and antiviral efficacy. Alternatively, the manual cleaning can serve as preparation for mechanical cleaning.

Tools:

- ALPRO® WL-clean cleaning solution
- ALPRO® WL-cid disinfectant solution
- ALPRO® WL-dry drying spray
- ALPRO® cleaning adapter, depending on the connector
- Lint-free disposable cloth

8.2.1 Manual cleaning

Connect the adapter to the Alpro® WL-clean spray bottle and insert it into the handpiece on the motor side. Clean the instrument with 3 spray puffs for 2 seconds and with 1 minute residence time.

8.2.2 Manual disinfection, rinsing and drying

Disinfect the outer surfaces of the instrument using a disposable cloth with Alpro® WL-cid. Connect the adapter to the Alpro® WL-cid spray bottle, insert it into the handpiece on the motor side, then spray for 3 seconds. After 5 minutes of residence time, dry the instrument with Alpro® WL-dry spray.

CAUTION

The residence time after the last spray puff is 5 minutes.

8.3 Automatic processing

8.3.1 Automatic cleaning, disinfection and rinsing

NOTE

NOTE

MK-dent recommends thermal disinfection devices that satisfy the requirements of EN ISO 15883-1 and are approved by the manufacturer for dental handpieces and turbines. Cleaning should be performed at min. 55 °C (131 °F) for at least 6 minutes and disinfection at min. 90 °C (194 °F) for at least 5 minutes (for A0 value > 3000). A mildly alkaline cleaning agent with pH 9 to 11, such as Dr. Weigert Neodisher® MediClean forte is recommended.

CAUTION

For automatic cleaning of the interior, use the adapters specified by the device manufacturer.

Only process instruments that are dry and free of residue using an automatic cleaning and disinfection process.

The automatic cleaning and disinfection process was validated using the device WD BHT INNOVA® M3 in program no. 03 (Thermal disinfection), with Dr. Weigert Neodisher® MediClean forte 0.5% as the cleaning agent. The instruments were connected with a silicon hose with a Luer-Lock adapter to the machine.

Parameters for automatic cleaning and disinfection:

- Pre-rinse for 1 minute with Neodisher® MediClean forte 0.5% cleaning agent at 30 °C (86 °F)
- Clean for 6 minutes with Neodisher® MediClean forte 0.5% cleaning agent at 55 °C (131 °F)
- Rinse for 1 minute with demineralized water
- Carry out thermal disinfection for 5 minutes at 90 °C (194 °F) using low-germ (< 10 CfU/ml) and low-endotoxin, (< 0.25 endotoxin units/ml) demineralized water (for an A0 value > 3,000).

8.3.2 Automatic drying

The drying process is normally part of the cleaning process of your thermal disinfectant. Please follow the instructions for your thermal disinfectant. Allow the instrument to cool down before removing it from the adapter to ensure this is safe and easier to do.

If there is still visible contamination in the thermal disinfectant, the process must be repeated.

To avoid any kind of damage to the instrument, make sure after each cycle that the instrument is dry inside and outside and free of residue. If moisture remains on or in the instrument after the cleaning process, dry it with a lint-free cloth. Make sure that the cavities of the instruments are dry and dry them with medical compressed air if necessary.

8.3.3 Maintenance, checks and testing

CAUTION

Premature wear and malfunctions from improper maintenance and care
 → Reduced product service life
 → Perform appropriate maintenance and care at regular intervals


NOTE

NOTE

The instrument requires maintenance after each processing and before each sterilization. Use MK-dent Premium Service Oil LU1011 for this purpose.

The instrument must be maintained after each processing and before each sterilization. Carry out the following steps in the specified order:

- Screw the included LT1014P oil adapter for maintenance of the spindle onto the oil can. Spray into the spindle for 2 seconds.
- Screw the included LT1013P oil adapter for the motor connection onto the oil can. Insert it into the instrument on the motor side and spray into the connection for 2 seconds.

8.4 Sterilization in an autoclave

8.4.1 Packaging

NOTE

NOTE

Use a validated packaging procedure that is suitable for sterilization. The sterilization pouch must be big enough for the instrument so that the packaging does not get stretched. The sterilization packaging must comply with EN ISO 11607-1 and be suitable for the sterilization process. Pack each instrument individually in a simple sterilization pouch. Seal the sterilization pouch carefully and check the seal for leaks.

8.4.2 Sterilization in a steam sterilizer according to EN 13060 / EN ISO 17665

WARNING

Premature wear and malfunctions from improper maintenance and care
 → Reduced product service life
 → Maintenance of the instrument with MK-dent Premium Service Oil is required before each sterilization.

Sterilization must be carried out before the first use and after each patient to prevent direct and cross-infection.

To avoid infection, wear suitable protective clothing such as a mask, gloves and safety goggles.

Contact corrosion from moisture can cause damage to the product.

→ Take the instrument out of the steam sterilizer directly after the sterilization cycle.

This medical device can withstand temperatures up to max. 135 °C (275 °F).

CAUTION

Depending on the conditions of use of your sterilizer or your sterilization methods, this product may undergo discoloration.

NOTE

NOTE

The sterilizer manufacturer's instructions for operation and load configuration must be followed explicitly. When sterilizing multiple instruments in one sterilization cycle, make sure that the instruments do not touch each other to ensure the best possible steam penetration. Observe the information in the sterilizer manufacturer's instructions for use.

After sealing it in the sterilization pouch, put the instrument onto a sterilization tray to avoid direct contact with the inner chamber surface of the sterilizer and then place it into the autoclave. MK-dent recommends the use of a sterilizer in accordance with EN 13060, EN ISO 17665 or EN 285.

8.4.3 Sterilization parameters

Only use the following sterilization parameters to achieve an SAL of 10⁻⁶:

Autoclave with a triple pre-vacuum (recommended): at least 5 minutes at min. 134 °C (273 °F).

The sterilization was validated with the device Tuttnauer 3870 EHS, type B autoclave, in a double pre-vacuum process at 134 °C (273 °F) for 2.5 minutes (half-cycle) with a drying time of 10 minutes. The instruments have been double-bagged in sterilization pouches and sealed.

The instructions listed above for processing a medical device for reuse have been validated and found to be suitable by the medical device manufacturer. The user is responsible for ensuring that the processing actually performed, including the devices and materials used in the processing facility as well as the personnel involved, will achieve the intended result. This requires verification and/or validation as well as routine monitoring of the process.

8.4.4 After sterilization

When the sterilization cycle is finished, take the instrument out of the autoclave. When doing so, be sure to use appropriate personal protective equipment, as the instruments will be hot.

8.5 Storage

After sterilization and drying, check that there are no drops of residual water. Store the processed instruments in a dry, dark and cool room, protected from dust and with the lowest possible microbial load.

Temperature	-20 °C to +70 °C (-4 °F to +158 °F)
Relative humidity	5 % - 95 %, non-condensing
Air pressure	700 hPa - 1060 hPa (10 psi - 15 psi)

CAUTION

Instrument malfunction during operation after heavily cooled storage conditions
→ Allow heavily cooled instruments to warm up to room temperature before starting to use them.

The instrument must be completely dried before storage. Any remaining water can result in discoloration or rust.

Do not store the instrument in a place where it could be exposed to water.

Do not store the instrument in places where chemicals are stored or where corrosive gases exist.

9 Maintenance

9.1 Regular maintenance and checks

WARNING

Inappropriate maintenance or service can result in malfunctions and/or deterioration.

→ Always conduct proper maintenance and servicing according to section **8.3.3 Maintenance, checks and testing**

CAUTION

For safe use of the product, perform checks and maintenance of the instrument according to these instructions for use.

9.2 Cleaning in case of problems

In case of problems with the air or water supply lines, clean each part of the instrument. The instrument must be disconnected from the hose coupling for cleaning. Instructions can be found in section **8.1 Pre-treatment at the point of use**.

10 Transport and disposal

CAUTION

Users or third parties may be infected by contaminated medical devices. Process the product before transport or disposal or dispose of it properly in compliance with the applicable regulations for the disposal of contaminated medical waste.

Process the instrument in accordance with section **8 Processing according to ISO 17664-1** and return it to us in the original packaging. We will take care of the professional disposal of the device and the original packaging for you.

11 Troubleshooting

NOTE

NOTE

This instrument must not be repaired by the user in the event of a malfunction. Please contact your dealer or the MK-dent repair service.

If this instrument does not function as described in these instructions for use, check the list below for your problem. If your problem is not resolved or if your problem is not listed below, please contact your dealer.

Symptom	Possible causes	Countermeasures
1) No water	- The water supply system of your treatment unit is switched off - No water gets into the hose - Water tubes are clogged - O-ring(s) on your motor is (are) damaged	- Activate the water supply system of your treatment unit - Check the water supply system of your treatment unit - Check the hose for potential contamination - Replace the O-ring(s) on your motor - Request repair by your dealer
2) Water leak on the hose connection	O-ring(s) on your hose connection is (are) damaged.	- Replace the O-ring(s) on your hose connection. - Request repair by your dealer
3) No light	- Glass rod is defective - Contactor of the light guide is defective	- Replace the tube - Repair/replace the connector - Request repair by your dealer
4) Irregular operating noises	- Ball bearing is dry, not lubricated sufficiently - Ball bearing defective - Corrugated spring washer is worn out - O-rings are worn out	- Lubricate/replace angle handpiece - Replace the corrugated spring washer - Replace O-rings - Request repair by your dealer
5) Strong vibrations	- Ball bearing defective - Ball bearing worn out	Request repair by your dealer.
6) Strong heat development	- Handpiece is dry, not lubricated sufficiently - O-rings are worn out - Ball bearing defective	- Lubricate the instrument - Replace O-rings - Request repair by your dealer
7) Optical damages	- Drop damage - Deformation	- Replacement - Request repair by your dealer

Symptom	Possible causes	Countermeasures
8) Rotation speed inadequate	• Air supply to the air-powered motor is insufficient • Supply voltage of the e-motor is too low • Ball bearing defective	• Check the compressed air supply • Check the supply voltage • Request repair by your dealer

12 Repair

CAUTION

Contaminated medical devices can infect users or third parties. Process the instrument before dispatching it.

Process this instrument in accordance with section **8 Processing according to ISO 17664-1** and send it to your dealer or directly to MK-dent if it is defective. Please contact MK-dent customer service to arrange for shipping or collection for factory repair. Questions about noise development, maintenance, preparation and functionality can often be clarified in a conversation, making it unnecessary to send in the product.

13 Technical data

Item no.	Type	Transmission/ reduction	Marking	Max. rotation speed [min ⁻¹]	Max. motor speed [min ⁻¹]	Coating	Sound level [dB(A)]
LB01	Straight handpiece	1:1	blue	20,000	20,000	Chrome	≤ 70
LB01W	Straight handpiece	1:1	blue	20,000	20,000	Chrome	≤ 70
LB02	Straight handpiece	1:1	blue	20,000	20,000	Chrome	≤ 70
LB1111	Contra-angle handpiece	1:1	blue	20,000	20,000	Chrome	≤ 70
LB1111C	Contra-angle handpiece	1:1	blue	20,000	20,000	Chrome	≤ 70
LB11W	Contra-angle handpiece	1:1	blue	20,000	20,000	Chrome	≤ 70
LB4111	Contra-angle handpiece	4:1	green	5,000	20,000	Chrome	≤ 70
LB4111P	Contra-angle prophyl handpiece	4:1	green	5,000	20,000	Chrome	≤ 70
LBH11	Contra-angle head	1:1	blue	20,000	20,000	Chrome	-
LBH11C	Contra-angle head	1:1	blue	20,000	20,000	Chrome	-
LBH11P	Contra-angle prophyl head	1:1	blue	20,000	20,000	Chrome	-
LE01	Straight handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LE01L	Straight handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LE11	Contra-angle handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LE11L	Contra-angle handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LE148	Contra-angle handpiece	14.8:1	green	2,703	40,000	Chrome	≤ 70
LE148L	Contra-angle handpiece	14.8:1	green	2,703	40,000	Chrome	≤ 70
LE148LP	Contra-angle prophyl handpiece	14.8:1	green	2,703	40,000	Chrome	≤ 70
LE148P	Contra-angle prophyl handpiece	14.8:1	green	2,703	40,000	Chrome	≤ 70
LE15	Contra-angle handpiece	1:5	red	200,000	40,000	Chrome	≤ 70
LE15L	Contra-angle handpiece	1:5	red	200,000	40,000	Chrome	≤ 70
LE16	Contra-angle handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LE16L	Contra-angle handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LE21	Contra-angle handpiece	2:1	green	20,000	40,000	Chrome	≤ 70
LE21L	Contra-angle handpiece	2:1	green	20,000	40,000	Chrome	≤ 70
LE21P	Contra-angle prophyl handpiece	2:1	green	20,000	40,000	Chrome	≤ 70
LE27	Contra-angle handpiece	2.7:1	green	14,814	40,000	Chrome	≤ 70
LE27L	Contra-angle handpiece	2.7:1	green	14,814	40,000	Chrome	≤ 70
LE54	Contra-angle handpiece	5.4:1	green	7,407	40,000	Chrome	≤ 70
LE54L	Contra-angle handpiece	5.4:1	green	7,407	40,000	Chrome	≤ 70
LE54LP	Contra-angle prophyl handpiece	5.4:1	green	7,407	40,000	Chrome	≤ 70
LE54P	Contra-angle prophyl handpiece	5.4:1	green	7,407	40,000	Chrome	≤ 70
LE74	Contra-angle handpiece	7.4:1	green	5,406	40,000	Chrome	≤ 70
LE74L	Contra-angle handpiece	7.4:1	green	5,406	40,000	Chrome	≤ 70
LC01	Straight handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LC01L	Straight handpiece	1:1	blue	40,000	40,000	Chrome	≤ 70
LC11	Contra-angle body	1:1	blue	40,000	40,000	Chrome	-
LC11L	Contra-angle body	1:1	blue	40,000	40,000	Chrome	-
LC27	Contra-angle body	2.7:1	green	14,814	40,000	Chrome	-
LC27L	Contra-angle body	2.7:1	green	14,814	40,000	Chrome	-
LC74	Contra-angle body	7.4:1	green	5,406	40,000	Chrome	-
LC74L	Contra-angle body	7.4:1	green	5,406	40,000	Chrome	-
LCH11	Contra-angle head	1:1	blue	40,000	40,000	Chrome	-
LCH16	Contra-angle head	1:1	blue	40,000	40,000	Chrome	-
LCH21	Contra-angle head	2:1	green	20,000	40,000	Chrome	-
LCH21P	Contra-angle prophyl head	2:1	green	20,000	40,000	Chrome	-
LP01L	Straight handpiece	1:1	blue	40,000	40,000	Titanium	≤ 70

Item no.	Type	Transmission/ reduction	Marking	Max. rotation speed [min ⁻¹]	Max. motor speed [min ⁻¹]	Coating	Sound level [dB(A)]
LP11L	Contra-angle handpiece	1:1	blue	40,000	40,000	Titanium	≤ 70
LP15L	Contra-angle handpiece	1:5	red	200,000	40,000	Titanium	≤ 70
LP27L	Contra-angle handpiece	2.7:1	green	14,814	40,000	Titanium	≤ 70
LP54L	Contra-angle handpiece	5.4:1	green	7,407	40,000	Titanium	≤ 70
LP74L	Contra-angle handpiece	7.4:1	green	5,406	40,000	Titanium	≤ 70
LPH11	Contra-angle head	1:1	blue	40,000	40,000	Titanium	-
LPH16	Contra-angle head	1:1	blue	40,000	40,000	Titanium	-
LPH21	Contra-angle head	2:1	green	20,000	40,000	Titanium	-

14 Explanation of symbols

	Caution – see section 1 Explanation of the warnings
	Sterilizable in a steam sterilizer (autoclave) up to the specified temperature
	Suitable for automatic processing in thermal disinfection devices
	Consult instructions for use and refer to the source of the electronic instructions for use
	CE marking with identification number of the notified body – indicates that the product is in compliance with European legislation for medical devices
	Medical device (CE conformity assessment procedure)
	Catalog number
	Serial number
	Manufacturer
	Date of manufacture
	Unique product identification – UDI (Unique Device Identifier) – indicates a barcode as containing Unique Device Identification
	HIBC Code (Health Industry Bar Code)

15 Guarantee and warranty

MK-dent assures you that your instrument has undergone a thorough **quality inspection** before leaving the factory. These quality inspections are carried out according to considerably stricter standards than the requirements of the product standard (ISO 14457 Dentistry - Handpieces and motors). In addition, we offer you comprehensive warranty services to ensure your satisfaction.

15.1 Guarantee

MK-dent grants the end customer a guarantee for perfect function and freedom from defects in material and workmanship from the date of purchase.

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (Instruments must be packed individually and with sufficient padding)
- No enclosure of the invoice copy of the device purchase (invoice recipient must be a dental practice)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

Regions, time periods and validity:

International	
Basic Line	12 months
Eco Line	24 months
Classic Line	24 months
Prime Line	24 months

15.2 Legal warranty

MK-dent grants the statutory warranty on all products to the contractual partner. This covers material and processing defects that were already present at the time of the transfer of risk or that occur within the statutory warranty period. Warranty claims will only be recognized if the product is submitted with proof of purchase in the form of an invoice or a copy of the delivery note (contractual partner, date of purchase, model and serial number must be clearly visible).

Exclusions:

- Wear parts, e.g., collets, ball bearings, gaskets, adapters, bulbs and O-rings
- Fall damage
- Willful damage and improper use contrary to the instructions for use
- Transport damage due to improper packaging (instruments must be packed individually and with sufficient padding)
- The warranty is valid from the date of purchase and is non-transferable
- Defects or their consequences due to interventions or modifications to the product by the end user or by a third party not authorized by MK-dent

15.3 Conditions

- Proof of purchase (invoice or delivery note) is required for guarantee or warranty claims.
- The guarantee or warranty is void if the product is tampered with or modified by unauthorized third parties.
- Repairs or replacements under warranty do not entitle the warranty period to be extended or restarted.
- MK-dent reserves the right to repair or replace defective instruments.

15.4 Utilization

In the event of a warranty claim, please contact MK-dent customer service directly. In the event of a warranty claim, please contact your contractual partner.

Contact MK-dent customer service

Phone: +49 4532 40049-0

Email: service@mk-dent.com

Please have the following information ready:

- Product model and serial number
- Proof of purchase with date
- Detailed description of the problem

Our customer service team will guide you through the warranty process and ensure a quick and satisfactory solution. After the repair, your instrument will be returned to you securely packaged and sealed together with a detailed test report.

15.5 Exclusion of liability

Within the scope of the guarantee or warranty, MK-dent is not liable for financial losses, downtimes, loan or rental equipment, loss of profit or similar consequences. The liability of MK-dent is limited to the purchase value of the product.

MK-dent is not liable for any defects and their consequences that may arise due to natural wear and tear, improper handling, improper cleaning or maintenance, non-compliance with operating or handling instructions, calcination or corrosion, contaminated air and water supply or chemical or electrical influences that are unusual or not permitted according to the instructions for use of MK-dent or other manufacturer's instructions for use.

Any liability is excluded if defects or their consequences are caused by interventions or modifications to the product by the end user or by a third party not authorized by MK-dent.

MK-dent reserves the right to change the warranty conditions at any time. Changes have no retroactive validity for products already purchased.

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
夹持力测量工具 RT1020	 A black handheld digital force gauge with a small LCD screen and four buttons labeled 'UNIT', 'TARE', 'LIGHT', and 'ON/OFF'. The brand name 'MK-dent' is printed on the device.
器械支架 RT2050	 A grey plastic base with multiple vertical slots holding several cylindrical metal handpieces of different sizes. The brand name 'MK-dent' is printed on the base.
MK-dent Premium Service Oil LU1011	 A blue aerosol can of premium service oil. The label features the 'MK-DENT' logo, a drop icon, and the text 'PREMIUM SERVICE OIL'. Below this, it says 'LU1011 650' and 'Volume = 500ml (16.9 fl. oz.)'. There are also three small icons of handpieces at the bottom of the label.

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
5.1 Acessórios

Ferramenta de medição da força de retenção RT1020	 A black handheld device with a digital display and four buttons labeled 'UNIT', 'TARE', 'LIGHT', and 'ON/OFF'. The MK-dent logo is visible at the bottom of the device.
Suporte de instrumentos RT2050	 A grey plastic holder with multiple slots for holding dental instruments. The MK-dent logo is visible on the front.
Premium Service Oil LU1011 da MK-dent	 A blue aerosol can with a white cap. The label features the MK-DENT logo, a drop icon, and the text 'PREMIUM SERVICE OIL'. Below this, it says 'LU1011' and '650'. At the bottom, it indicates 'Volumen = 500ml (l)' and 'NET 12.21 oz (345g)'. There are also three small icons of dental handpieces.

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